

TOPSAFE HEADREST COVER

HEADREST COVER MADE OF STERILIZATION PAPER



Quality:	crepe, ca. 58 g/m ² +/- 2 wood free, glued, made of pure cellulose
Colors:	white, green, blue, yellow
Sizes:	280 x 360 mm 200 pcs/carton
Packaging control:	electronical control system for sheet-quantity per carton, stacked, tied up and wrapped in stretch foil on palettes
Manufacturing control:	according to DIN EN 868-2; DIN EN 868-1 and DIN 58 953, 2.14 + 2.15
Shelf life:	minimum six months unlimited under following storage conditions:
Storage conditions:	no extremely fluctuation of temperature in the warehouse, temperature between 10 and 30 °C No storage in direct sun, dry storage