


SILK



	<i>Silk Braided</i>	<i>Silk Virgin</i>
Type	Braided multifilament	Twisted multifilament
Composition	Braided fibres from the cocoon of the silkworm	Twisted fibres from the cocoon of the silkworm
Coating	Waxes Silicone	None
Colour	Black	Blue
Tissue reaction	Moderate	
Absorption	Silk suture elicits an initial inflammatory reaction in tissues, which is followed by gradual encapsulation of the suture by fibrous connective tissues.	
Presentation	Box with 12 sutures	

*For quality
&
safety*

smi

www.sutures.be

SILK

CHARACTERISTICS

- ✓ High tensile strength, non absorbable
- ✓ Braided or twisted structure

- ✓ Coated multifilament

- ✓ Hermetically sealed packing

BENEFITS

- Good and extended tissue support for up to 3 month

- Excellent handling properties
High flexibility
High tensile strength
Excellent knot security

- Soft passage through the tissues with minimal sawing, tissue drag and trauma
Good knot tie down / adjustability
Reduced capillary action

- Guaranteed seal and product sterility

INDICATIONS General soft tissue closing and/or ligation
especially in General surgery
Skin closure
Ophthalmic surgery
Oral surgery
Gastrointestinal surgery
Gynaecology
Obstetrics

CONTRAINDICATIONS Urinary tract
Biliary tract

Known allergies or sensitivities to silk

Without needle



Metric **USP** **Length** **Order number**

2	3/0	150 cm	<i>black</i>	8200015
3	2/0	150 cm	<i>black</i>	8300015
3.5	0	150 cm	<i>black</i>	8350015
4	1	150 cm	<i>black</i>	8400015
5	2	150 cm	<i>black</i>	8500015
6	3&4	150 cm	<i>black</i>	8600015
8	6	150 cm	<i>black</i>	8800015
2	3/0	6 x 45 cm	<i>black</i>	8200044
3	2/0	6 x 45 cm	<i>black</i>	8300044
3.5	0	6 x 45 cm	<i>black</i>	8350044
4	1	6 x 45 cm	<i>black</i>	8400044
5	2	6 x 45 cm	<i>black</i>	8500044
1.5	4/0	17 x 45 cm	<i>black</i>	8150046
2	3/0	17 x 45 cm	<i>black</i>	8200046
3	2/0	17 x 45 cm	<i>black</i>	8300046
3.5	0	17 x 45 cm	<i>black</i>	8350046
4	1	17 x 45 cm	<i>black</i>	8400046
5	2	17 x 45 cm	<i>black</i>	8500046
1.5	4/0	13 x 60 cm	<i>black</i>	8150060
2	3/0	13 x 60 cm	<i>black</i>	8200060
3	2/0	13 x 60 cm	<i>black</i>	8300060
3.5	0	13 x 60 cm	<i>black</i>	8350060
4	1	13 x 60 cm	<i>black</i>	8400060
2	3/0	12 x 75 cm	<i>black</i>	8200075
3	2/0	12 x 75 cm	<i>black</i>	8300075
3.5	0	12 x 75 cm	<i>black</i>	8350075
4	1	12 x 75 cm	<i>black</i>	8400075
5	2	12 x 75 cm	<i>black</i>	8500075

Without needle



Ligapack



Metric **USP** **Length** **Order number**

2	3/0	250 cm	<i>black</i>	8200025
3	2/0	250 cm	<i>black</i>	8300025
3.5	0	250 cm	<i>black</i>	8350025
4	1	250 cm	<i>black</i>	8400025
5	2	250 cm	<i>black</i>	8500025
6	3&4	250 cm	<i>black</i>	8600025

Round bodied taper point



1/2 circle



		Metric	USP	Length	Color	Order number
HR 12 mm		0.7	6/0	45 cm	black	8070112
		1	5/0	45 cm	black	8100112
		1.5	4/0	75 cm	black	8150112
HR 16 mm		0.7	6/0	45 cm	black	8070116
		1	5/0	75 cm	black	8100116
		1.5	4/0	75 cm	black	8150116
HR 20 mm		2	3/0	75 cm	black	8200116
		1	5/0	75 cm	black	8100120
		1.5	4/0	75 cm	black	8150120
HR 26 mm		2	3/0	75 cm	black	8200120
		3	2/0	75 cm	black	8300120
		3.5	0	75 cm	black	8350120
		1	5/0	75 cm	black	8100126
		1.5	4/0	75 cm	black	8150126
HR 30 mm		2	3/0	75 cm	black	8200126
		3	2/0	75 cm	black	8300126
		3.5	0	75 cm	black	8350126
		4	1	75 cm	black	8400126
		5	2	75 cm	black	8500126
		2	3/0	75 cm	black	8200130
HR 37 mm		3	2/0	75 cm	black	8300130
		3.5	0	75 cm	black	8350130
		4	1	75 cm	black	8400130
		5	2	75 cm	black	8500130
HR 40 mm		2	3/0	75 cm	black	8200137
		3	2/0	75 cm	black	8300137
		3.5	0	75 cm	black	8350137
HR 45 mm		4	1	75 cm	black	8400137
		5	2	75 cm	black	8500137
		3.5	0	75 cm	black	8350140
HR 45 mm		4	1	75 cm	black	8400140
		5	2	75 cm	black	8500140
		3.5	0	75 cm	black	8350145
HR 45 mm		4	1	75 cm	black	8400145
		5	2	75 cm	black	8500145
		3.5	0	75 cm	black	8350145

Round bodied taper point



		Metric	USP	Length	Order number
--	--	--------	-----	--------	--------------

3/8 circle



DR 12 mm		1	5/0	45 cm	<i>black</i> 8100212
DR 16 mm		1.5	4/0	75 cm	<i>black</i> 8150216
		2	3/0	75 cm	<i>black</i> 8200216
DR 25 mm		2	3/0	75 cm	<i>black</i> 8200225
		3	2/0	75 cm	<i>black</i> 8300225
DR 30 mm		3	2/0	75 cm	<i>black</i> 8300230

Round bodied taper point



		Metric	USP	Length	Order number
--	--	--------	-----	--------	--------------

Straight

GR 40 mm		3	2/0	75 cm	<i>black</i> 8302040
----------	--	---	-----	-------	----------------------

Reverse cutting



		Metric	USP	Length	Order number
--	--	--------	-----	--------	--------------

3/8 circle



DS 12 mm		0.5	7/0	45 cm	<i>black</i> 8051512
		0.7	6/0	45 cm	<i>black</i> 8071512
		1	5/0	45 cm	<i>black</i> 8101512
		1.5	4/0	75 cm	<i>black</i> 8151512
		2	3/0	75 cm	<i>black</i> 8201512
DS 12 mm Easy Pass		0.7	6/0	50 cm	<i>black</i> 8071512P
DS 14 mm		0.7	6/0	45 cm	<i>black</i> 8071514
DS 16 mm		0.7	6/0	75 cm	<i>black</i> 8071516
		1	5/0	75 cm	<i>black</i> 8101516
		1.5	4/0	75 cm	<i>black</i> 8151516
		2	3/0	75 cm	<i>black</i> 8201516
		3	2/0	75 cm	<i>black</i> 8301516
DS 16 mm Easy Pass		1	5/0	50 cm	<i>black</i> 8101516P
		1.5	4/0	50 cm	<i>black</i> 8151516P
DS 16 mm black		1	5/0	75 cm	<i>black</i> 8101516VB
		1.5	4/0	75 cm	<i>black</i> 8151516VB

Reverse cutting



Metric

USP

Length

Order number

DS 19 mm



0.7	6/0	75 cm	black	8071519
1	5/0	75 cm	black	8101519
1.5	4/0	75 cm	black	8151519
2	3/0	45 cm	black	8201519-45
2	3/0	75 cm	black	8201519
3	2/0	75 cm	black	8301519
3.5	0	75 cm	black	8351519
4	1	75 cm	black	8401519

DS 19 mm
Easy Pass



1.5	4/0	50 cm	black	8151519P
2	3/0	50 cm	black	8201519P

DS 19 mm black



1	5/0	75 cm	black	8101519VB
1.5	4/0	75 cm	black	8151519VB
2	3/0	75 cm	black	8201519VB

DS 24 mm



1	5/0	75 cm	black	8101524
1.5	4/0	75 cm	black	8151524
2	3/0	75 cm	black	8201524
3	2/0	75 cm	black	8301524
3.5	0	75 cm	black	8351524
4	1	75 cm	black	8401524

3/8 circle



DS 24 mm
Easy Pass



2	3/0	50 cm	black	8201524P
---	-----	-------	-------	----------

DS 26 mm



2	3/0	75 cm	black	8201526
3	2/0	75 cm	black	8301526

DS 30 mm



1.5	4/0	75 cm	black	8151530
2	3/0	75 cm	black	8201530
3	2/0	75 cm	black	8301530
3.5	0	75 cm	black	8351530
4	1	75 cm	black	8401530
5	2	75 cm	black	8501530

DS 37 mm







2	3/0	75 cm	black	8201537
3	2/0	75 cm	black	8301537
3.5	0	75 cm	black	8351537
4	1	75 cm	black	8401537
5	2	75 cm	black	8501537
6	3&4	75 cm	black	8601537


DS 45 mm











3	2/0	75 cm	black	8301545
3.5	0	75 cm	black	8351545
4	1	75 cm	black	8401545
5	2	75 cm	black	8501545




Reverse cutting  **Metric** **USP** **Length** **Order number**





3/8 circle 	DS 60 mm 	4	1	75 cm	<i>black</i>	8401560
	DS 90 mm 	3.5	0	100 cm	<i>black</i>	8351590




Reverse cutting  **Metric** **USP** **Length** **Order number**

1/2 circle 	HS 12 mm 	0.7	6/0	75 cm	<i>black</i>	8071412
	HS 17 mm 	1	5/0	75 cm	<i>black</i>	8101417
		1.5	4/0	75 cm	<i>black</i>	8151417
		2	3/0	75 cm	<i>black</i>	8201417
	HS 20 mm 	1.5	4/0	75 cm	<i>black</i>	8151420
		2	3/0	75 cm	<i>black</i>	8201420
		1	5/0	75 cm	<i>black</i>	8101422
	HS 22 mm 	1.5	4/0	75 cm	<i>black</i>	8151422
		2	3/0	75 cm	<i>black</i>	8201422
		3	2/0	75 cm	<i>black</i>	8301422
	HS 26 mm 	1	5/0	75 cm	<i>black</i>	8101426
		1.5	4/0	75 cm	<i>black</i>	8151426
		2	3/0	75 cm	<i>black</i>	8201426
		3	2/0	75 cm	<i>black</i>	8301426
	HS 36 mm 	3.5	0	75 cm	<i>black</i>	8351426
3.5		0	75 cm	<i>black</i>	8351436	
4		1	75 cm	<i>black</i>	8401436	

Cutting  **Metric** **USP** **Length** **Order number**

Straight 	GS 40 mm 	3	2/0	75 cm	<i>black</i>	8302140
	GS 60 mm 	2	3/0	75 cm	<i>black</i>	8202160
		3	2/0	75 cm	<i>black</i>	8302160
		3.5	0	75 cm	<i>black</i>	8352160
		4	1	75 cm	<i>black</i>	8402160

Round bodied taper cutting point 		Metric	USP	Length	Order number	
1/2 circle 	HRT 18 mm 	1	5/0	75 cm	<i>black</i>	8100518
		1.5	4/0	75 cm	<i>black</i>	8150518
		2	3/0	75 cm	<i>black</i>	8200518
	HRT 27 mm 	3	2/0	75 cm	<i>black</i>	8300518
		2	3/0	75 cm	<i>black</i>	8200527
		3	2/0	75 cm	<i>black</i>	8300527

Round bodied taper cutting point 		Metric	USP	Length	Order number	
3/8 circle 	DRT 18 mm 	1.5	4/0	75 cm	<i>black</i>	8150418
		2	3/0	75 cm	<i>black</i>	8200418