

Two Striper®



Two Striper® Carbopost®

Two Striper® Carbopost® is a fiber reinforced composite post that's as strong as metal, which reduces shear stress.

Made of 60% T700 carbon fibers and 40% epoxy resin.

- Autoclavable up to 132°C/ 269.6°F
- Bio-compatible and corrosion free
- Anatomical shape and 5 different diameters (0.8mm – 1.6mm)



Mechanical Properties	
	Two Striper® Carbopost®
Elastic Modulus	141GPa
Shear Strength	23.9 MPa



Two Striper® Carbopost® Kits			
C801001	Basic Kit : 50 assorted posts (1.0, 1.2, 1.4, 1.6) and 4 matching finishing drills		
C801002	Intro Kit: 20 assorted posts (1.0, 1.2, 1.4, 1.6) and 4 matching finishing drills		
Two Striper® Carbopost® Refills			
10 Posts/Pack		10 posts/pack with matching finishing drill	
C802013	Ø 0.8		
C802001	Ø 1.0	C802006	Ø 1.0
C802002	Ø 1.2	C802007	Ø 1.2
C802003	Ø 1.4	C802008	Ø 1.4
C802004	Ø 1.6	C802009	Ø 1.6