

**VASO04-VASO RECOLECTOR AUTOCLAVABLE**

Vaso recolector de aspiración sin protector de rebose. Autoclavable.

VASO04-REUSABLE COLLECTOR JAR

Reusable collector aspiration jar without protector of overflow.

**CARACTERISTICAS TECNICAS/TECHNICAL CHARACTERISTICS:**

CLASE/CLASS:	IIA
NORMATIVA/STANDARD:	UNE-EN ISO 10079-3
CAPACIDAD/CAPACITY:	225 cc
INTERVALO RECOMENDADO DE USO/RECOMMENDED USE INTERVAL:	3-6 bar
DEPRESION MAXIMA DE UTILIZACION/MAXIMUM USE DEPRESSION:	-740 mbar
DIMENSIONES/ DIMENSIONS:	177 x 86 x 64 mm
PESO/ WEIGHT:	170 gr
CAUDAL MAXIMO DE USO/MAXIMUM USE FLOW:	20 l/min
CONEXIÓN DE ENTRADA/INLET CONNECTION:	9 mm
CONEXIÓN DE SALIDA/OUTLET CONNECTION:	9 mm
REUTILIZABLE Y AUTOCLAVABLE/REUSABLE AND AUTOCLAVABLE:	130°C
CONDICIONES AMBIENTALES DE UTILIZACIÓN/ ENVIRONMENTAL CONDITIONS:	0°C +40°C
CONDICIONES AMBIENTALES DE ALMACENAMIENTO/ENVIRONMENTAL STORAGE CONDITIONS:	-20°C +60°C