



→ IDENTIFICATION

perio

→ MORPHOLOGY



102° angled tipholder with conical PEEK tip
(working length 10 mm, Ø apical 0,8 mm,
Ø coronal 2 mm)

→ SURFACE

brightly polished medical grade stainless steel
for tipholder ICS and 100% PEEK for IC1

→ TREATMENT

efficient and gentle cleaning of implants
and restorations



→ ICS



→ IC1

implant cleaning tip IC1
can be easily tightened
without special tools
and is fully autoclavable