



→ **MECTRON**
MICROPIEZO

NEW MICROPIEZO

New compact and essential design, maximum ergonomics: reduced fatigue for the operator



→ SOFT MODE

The ultra-gentle scaling for sensitive patients.

- mectron innovative SOFT MODE, allows a nearly painless treatment for patients and maximum comfort for clinicians.
- Available in all the applications

→ WiFi



Integrated WiFi technology - already prepared for future IoT functions.



→ **independent irrigation**

Independent bottle tank switching in different colour according to the clinical indication (capacity 500 ml).

→ **touch keyboard**

Exclusive multifunction touch keyboard with colour code according to the clinical indication.

→ **flush function**

Short flush cycle - to get the irrigation at the tip before starting the treatment; Long flush cycle - to perform an automatic cleaning cycle of the irrigation circuit at the end of each treatment.

→ **slim handpiece**

Slim ultrasonic LED handpiece with circular light – lightweight only 55 g!

→ **multifunctional**

3 functions with 12 different power settings
- scaling and perio
- endodontic
- prosthetics

→ **wide insert range**

- More than 30 inserts available
- 1 key 1 insert! Every insert is supplied with its own torque wrench to avoid cross contamination and to reduce the sterilization time.



SOFT MODE - THE ULTRA-GENTLE SCALING



“SOFT MODE, in this era of biofilm-focused care, is a critical tool to use at the beginning of preventive maintenance visits. I would recommend to use on adolescents and adults for professional biofilm removal to promote oral health and overall health of patients.”

Prof. Susan Wingrove, USA



“Mectron SOFT MODE allows the highest ultrasound efficiency avoiding the disadvantage of the pain for the patient ... it's a really comfortable treatment.”

Dr. Marie-Ève Bezzina, France



“Can be used in all the clinical situations, especially in case of particularly sensitive patients offering them a mininvasive treatment.”

Prof. Annamaria Genovesi, Italy



“It is more pleasant for patients and I like work with it because it is successful. It damaging far less of the dental tissue.”

Dr. Susanne Graack, Germany

mectron s.p.a., via Loreto 15/A,
16042 Carasco (Ge), Italy,
tel +39 0185 35361, fax +39 0185 351374
www.mectron.com, mectron@mectron.com



© Copyright mectron S.p.A., Carasco, Italy
All rights reserved. Texts, pictures and graphics of mectron brochures are protected by copyright and other protection laws. Without written approval of mectron S.p.A. the contents may not be copied, distributed, changed or made available to third parties for commercial purposes.