


	FICHA TÉCNICA	CODIGO: TDS_009
	GUANTES LÁTEX ESTÉRILES SIN POLVO	FECHA: 10/Abr/17
	BESTDENT	VERSION: 01

GUANTES ESTÉRILES LÁTEX SIN POLVO BESTDENT		
REFERENCIAS	74506: Tamaño 6 74507: Tamaño 6.5 74508 : Tamaño 7 74509: Tamaño 7.5 74510: Tamaño 8 74511: Tamaño 8.5	
FABRICANTE / RESPONSABLE	Anhui Haojie Plastic & Rubber Products Co. Ltd. Daqing Rd., 170 233010 Bengbu China	
CLASIFICACIÓN	Producto Sanitario Clase IIa, Directiva 93/42/CEE	
DESCRIPCIÓN	Guantes estériles de látex sin polvo, de un solo uso. Ambidiestros.	
ESPECIFICACIONES	Tamaños: 6, 6.5, 7, 7.5, 8, 8.5	
EMBALAJE	Estuche con 50 unidades	
CADUCIDAD	5 años	
FINALIDAD	Protección de barrera de la mano del usuario para la realización de exámenes médicos o clínicos, procedimientos diagnósticos, terapéuticos y quirúrgicos, en odontología y también para fines de laboratorio.	
INFORMACIÓN SOBRE LA SEGURIDAD	Los guantes están fabricados con látex, que puede causar reacciones alérgicas en algunas personas. Si esto ocurre, discontinuar su uso inmediatamente y consultar con el médico.	
ALMACENAMIENTO Y TRANSPORTE	Almacenar en lugar fresco y seco. Mantener alejado de la luz solar. Proteger de la humedad y evitar el calor excesivo. Almacenar entre 10°C-35°C.	
NORMATIVAS QUE CUMPLEN	EN 455:-1,-2,-3,-4 EN 420+A1	

	TECHNICAL DATA SHEET	CODE: TDS_009
	STERILE LATEX GLOVES POWDER FREE	DATE: 10/Abr/17
	BESTDENT	VERSION: 01

STERILE LATEX GLOVES POWDER FREE BESTDENT		 bestdent GUANTES QUIRÚRGICOS ESTÉRILES SIN POLVO - 6 50 unidades <small> DE STERILE OP-HANDSCHUHE OHNE PUDER. EN STERILE SURGICAL GLOVES POWDER FREE. FR GANTS CHIRURGIE STÉRILES SANS POUFRE. IT GUANTI CHIRURGICI STERILI SENZA POLVERE. PT LUVAS CIRÚRGICAS ESTERILIZADAS SEM PÓ. </small>
REFERENCES	74506: size 6 74507: size 6.5 74508 : size 7 74509: size 7.5 74510: size 8 74511: size 8.5	
MANUFACTURER	Anhui Haojie Plastic & Rubber Products Co. Ltd. Daqing Rd., 170 233010 Bengbu China	
CLASSIFICATION	Medical Device, Class IIa, Directive 93/42/EEC	
DESCRIPTION	Sterile latex gloves powder free. Ambidextrous. For single use.	
SPECIFICACIONES	Sizes: 6, 6.5, 7, 7.5, 8, 8.5	
PACKAGING	50 units	
SHELF LIFE	5 years	
USE	Gloves intended for use in the medical field to protect patient and user from cross contamination. Intended for conducting medical examinations, diagnostic, therapeutic and surgical procedures, in dentistry and in handling contaminated medical material.	
SAFETY INFORMATION	The gloves are made of natural rubber latex that may cause allergic reactions in some people. If allergic reactions occur, discontinue use immediately and consult the doctor.	
HANDLING AND STORAGE CONDITIONS	Store in a cool dry place, away from direct sunlight. Protect from humidity and excessive heat. Storage between 10°C-35°C.	
STANDARDS APPLIED	EN 455:-1,-2,-3,-4 EN 420+A1	