


	TECHNICAL DATA SHEET	CODE: TDS_00X
	YELLOW ORAL TIPS	DATE: 14/March/20
	PROCLINIC	VERSION: 01

PRODUCT	YELLOW INTRA-ORAL TIPS PUNTAS INTRAORALES AMARILLAS	
REFERENCE	7680	
COMPANY	SHANGHAI NSJ Hardware	
CLASSIFICATION	CLASS I MEDICAL DEVICE Directive 93/42/CEE	
PRODUCT DESCRIPTION	YELLOW INTRA-ORAL TIPS PUNTAS INTRAORALES AMARILLAS	
SPECIFICATIONS	<p>Color: Yellow Connector ID: 2.1mm Connector OD: 3.6mm Outlet ID: 0.8mm Material: PP Matched Mixing Tip: 2# & 5#</p> <p>Color: amarillo ID del conector: 2.1 mm Conector OD: 3.6 mm ID de salida: 0.8mm Material: PP Punta de mezcla combinada: 2 # y 5 #</p>	
PACKAGING	50pcs/bag, 500bags/ctn 50 uds / bolsa, 500 bolsa / ctn	
SHELF LIFE	5 years 5 años	
INTENDED OF USE	Used together with 2# yellow & 5# blue mixing tips to precisely deliver the mixture of dental impression material to tiny specified area. Usually for the second impression.	

	TECHNICAL DATA SHEET	CODE: TDS_00X
	YELLOW ORAL TIPS	DATE: 14/March/20
	PROCLINIC	VERSION: 01

	Se utiliza junto con puntas de mezcla 2 # amarillas y 5 # azules para administrar con precisión la mezcla de material de impresión dental en un área específica pequeña. Por lo general, para la segunda impresión.
SECURITY	Yes Si
TRANSPORT AND STORAGE CONDITIONS	Normal temperature and pressure Temperatura y presión normal