



File No.:

Version: A

5. Endo Radar Pro 产品规格

Endo Radar Pro Product Specification

1、机器正常工作条件:

Working condition:

1.1、温度: $5^{\circ}\text{C} \sim 40^{\circ}\text{C}$ 。

1.1 Working temperature: $5^{\circ}\text{C} \sim 40^{\circ}\text{C}$ 。

1.2、相对湿度: 30%~75%。

1.2 Relative humidity: 30%~75%。

2、贮存条件:

产品应贮存在于相对湿度不超过 10%~93%,大气压力为 70KPa~106Kpa, $-20^{\circ}\text{C} \sim +55^{\circ}\text{C}$,通风良好的室内。

Storage conditions:

Relative humidity: 10%~93%

Atmosphere pressure: 70kPa~106kPa

Environment temperature: $-20^{\circ}\text{C} \sim +55^{\circ}\text{C}$

3、机器各参数:

Main technical specifications

3.1、电源: 输入: $\sim 100\text{V}-240\text{V}$ 50Hz/60Hz 0.4Amax

输出: DC5V/1A

Power supply: Input: $\sim 100\text{V}-240\text{V}$ 50Hz/60Hz 0.4Amax

Output:DC5V/1A



File No.:

Version: A

3.2、电池:

手柄马达锂电池: 3.7V /1200mAh

Battery:

Lithium Battery in motor handpeice: 3.7V /1200mAh

3.3、扭矩范围: 0.4Ncm-5.0Ncm (4mNm~50mNm)

Torque: 0.4Ncm-5.0Ncm (4mNm~50mNm)

3.4 转速范围: 100rpm~2500rpm

Rotate speed: 100rpm~2500rpm

3.5、按防电击类型分类: II 类

Type of protection against electric shock: Class II

3.6、按防电击程度分类: B 型应用部分

Degree of protection against electric shock: Type B applied part

3.7、按对进液的防护程度分类:普通器材(IPX0)

Degree of protection against harmful ingress of water: Ordinary equipment (IPX0)

3.8、易燃麻醉混合物和空气,氧气或一氧化二氮的安全度数:不适合在易燃麻醉混合物和空气,氧气或一氧化二氮中使用的器材。

Degree of safety of application in the presence of a flammable anesthetic being used in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.