
	TECHNICAL DATA SHEET	CODE: TDS_00X
	PINK MIXING TIPS	DATE: 18/May/18
	PROCLINIC	VERSION: 01

PRODUCT	PINK MIXING TIPS	
REFERENCE	78694	
COMPANY	SHANGHAI NSJ Hardware	
CLASSIFICATION	CLASS I MEDICAL DEVICE Directive 93/42/CEE	
PRODUCT DESCRIPTION	Pink mixing oral tips.	
SPECIFICATIONS	Color & Ratio: Pink 1:1 Total Length: 82mm Outlet Type and Size: 5.4mm Material: PP tube + POM element Matched Cartridge: 50ml 1:1 cartridge, 24ml 1:1 syringe Matched oral tip: N2, N5 intra-oral tip	
PACKAGING	50pcs/bag, 50bags/ctn	
SHELF LIFE	5 years	
INTENDED OF USE	For mixing various Regular/Heavy-body dental impression materials, direct hard relining materials, permanent soft A-silicone relining materials, bite registration materials etc.	
SECURITY	Yes	
TRANSPORT AND STORAGE CONDITIONS	Normal temperature and pressure	

