


DATASHEET FOR PHOTOGRAPHIC MIRRORS


Description : Double sided photographic mirrors

Dimensions :


MZP511 Ti
MZP512 Chr
MZP516 Pur



MZP541
MZP542
MZP546



MZP521
MZP522
MZP526




MZP551
MZP552
MZP556



MZP531
MZP532
MZP536




MZP581
MZP582
MZP586



MZP561
MZP562
MZP566



MZP571
MZP572
MZP576



GLASS:

Float glass
Thickness : 3mm
weight : 0,050 Kg

COATINGS :

CHROME COATING : through evaporation under vacuum
TITANIUM COATING : through evaporation under vacuum
PURE REFLECT COATING: Niobium/tantalum through evaporation under vacuum

CLASS : I

PACKAGING : plastic sleeve unit with protective foam