


Precision Machined GUTTA-PERCHA MASTER CONES



Precision machined
gutta-percha master cones

- Variable tapered design per
Tulsa Dental Specialties
file specification
- Predictably fits the
shaped root canal

Micronized formula

- Delivers warm gutta-percha
4.5+mm beyond the heat source
- Working Temperature of 105°-180° C

Latex Free
Made in the USA



1-800-662-1202

For the latest information consult
www.TulsaDentalSpecialties.com

DENTSPLY
TULSA DENTAL
SPECIALTIES