



TOP TRAY

DISPOSABLE MEDICAL TRAY PAPER

CONFIGURATION OVERVIEW

Material: Natural cellulose fibers
Type: Non-Sterile
 wood free
 Latex free



TECHNICAL SPECIFICATIONS

Optical brightener negative
Chlorine: negative












| Parameter | Standard | Data |
|--------------------------|------------------------|-------------------------|
| Basis weight | DIN EN ISO 536 | 55 ± 2 g/m ² |
| Thickness | DIN EN 20 534, ISO 534 | 0.16 ± 0.02 mm |
| Tensile strength MD | DIN EN ISO 1924-2 | ≥ 28 N/15 mm |
| Tensile strength MD, wet | DIN ISO 3781 | ≥ 4 N/15 mm |
| Stretch MD | DIN EN ISO 1924-2 | ≥ 4 % |

QUALITY SPECIFICATION

Conformity to standards: Produced according to ISO 9001:2015 and ISO 14001:2015

Storage: Store in a dry place at room temperature, away from heat sources, flames or sparks.

PRODUCT VARIATIONS

| | | | | | | |
|-------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| Color code | YE | LI | GR | PI | FU | LU |
| Color | Yellow | Lilac | Green | Pink | Fuchsia | Laguna Blue |
| |  |  |  |  |  |  |
| Color code | FG | BL | MB | WH | HO | |
| Color | Fresh Green | Blue | Med Blue | White | Orange | |
| |  |  |  |  |  | |

| Art-No. | Dimension |
|----------------------------|-----------|
| TRAAA04-NEU-Color Code-CA1 | 18x28 cm |
| TRAAA04-NEU-Color Code-C27 | 28x36 cm |

PACKAGING

Sheet/unit 250
Units/carton 10 (size 18x28 cm), 5 (size (28x36 cm)
Carton size 37.5x29.5x23 cm
Carton weight 7.3 kg
Pallets/carton 64
Carton/layer 8
Layers/pallet 8
Pallet size 80x120x198 cm
Pallet weight 487 kg