


Fi - P



0.2 second heating time

High-performance conversion system ensures the efficient cutting




0.35mm pre-bend tip

Thinner tips can penetrate deep into the narrow root canal, providing you a better treatment effect.



High precision temperature control

Constant temperature prevents scalding



4 working temperatures

150°C、180°C、200°C、230°C can match gutta-percha of different brands on the market.



4 hours of endurance

Large capacity battery & special heating system enable 1500 operations with a full charge.



Type	diameter	Taper
WP3504	0.35mm	0.04
WP4004	0.40mm	0.04
WP4504	0.45mm	0.04
WP5004	0.50mm	0.04
WP5506	0.55mm	0.06
WP5508	0.55mm	0.08
WP5510	0.55mm	0.10
WP6012	0.60mm	0.12
WP5508L	0.55mm	0.08

Fi-G:Time for 0°C-230°C (S)

Woodpecker

15

Korean B brand

24

Chinese C brand

51

Chinese D brand

56

Fi - G



Quick heating within 15 seconds

reduces time for waiting and increase your working efficiency.



0.5mm injection needle

Easy to reach the apical 3mm of canal brings better filling.



3mm



Prevent gutta-percha leakage


360 rotatable needle

Optimized structure, 360° uniform rotation, effectively avoiding gutta-percha leakage.



Pushing ram with anti-return structure

Mechanical-locking structure, effectively preventing the pushing ram from returning and avoiding the generation of cavities and bubbles.



Type	Size	Length
20G 22mm	20G	22mm
20G 24mm	20G	24mm
20G 28mm	20G	28mm
23G 24mm	23G	24mm
23G 28mm	23G	28mm
25G 24mm	25G	24mm