



Free Scan Imaging Plate Scanner



Guilin Woodpecker Medical Instrument Co.,Ltd.

Information Industrial Park,Guilin National High-Tech
Zone,Guilin,Guangxi,541004 P.R.China
Sales Dept: +86-773-5873196
E-mail: woodpecker@glwoodpecker.com
Website: <http://www.glwoodpecker.com>
ZMN-XC-074

Scan and log on to website
for more information



Find us on:




Find us on:



Find us on:



Bringing you a fully-controlled
and intelligent image viewing experience



Specifications

Imaging plate sizes: 0; 1; 2; 3

Grayscale: 65536 levels (16Bit)

Theoretical resolution: 17lp/mm

Connection mode: USB


Working mode: Online mode

Remarkable imaging sharpens each detail

25um scanning accuracy.

17LP/mm+65535 grayscale detail recognition ability.

Clear details and image processing make diagnosis and treatment more effective.



Easy to operate




One-button start for scanning.

120-second autolock
to keep data safe.

Advanced algorithm leads to intelligence

Multi-layered structural image algorithm framework
enhances the details of super-resolution images.

Multi-mode processing algorithm to meet the image viewing
habits of different doctors.






Quick image viewing

Imported ultra-thin imaging plate has moderate thickness to keep oral cavity comfortable.

The 0-3 size imaging plates are available for photographing all tooth positions.



Powerful chip Intelligent scanning

FPGA chip accelerates data acquisition.

Support image viewing on mobile phones, tablets, and computers.

Exclusive image viewing software.

Greatly facilitate image viewing and processing.

Real-time intelligent dose analysis Safe dental film photographing

Appropriate dose.

Accurately record the dose value for each photographing.


Timely reminder of exposure dose for clear images.

Ensure secure photographing.

Too low dose

Too high dose

Appropriate dose



Good-looking and handy Compacted size of 0.046m² and space-saving position

The silver-black scheme is suitable for upscale and official occasions.

The compacted size of 0.046 m² remains full-fledged Standing nearby armchair to safeguard your teeth

