


	TECHNICAL DATA SHEET	CODE: TDS_00X
	BROWN MIXING TIPS	DATE: 14/March/20
	PROCLINIC	VERSION: 01

PRODUCT	BROWN MIXING TIPS PUNTAS DE MEZCLA MARRÓN	
REFERENCE	97116	
COMPANY	SHANGHAI NSJ Hardware	
CLASSIFICATION	CLASS I MEDICAL DEVICE Directive 93/42/CEE	
PRODUCT DESCRIPTION	BROWN MIXING TIPS PUNTAS DE MEZCLA MARRÓN	
SPECIFICATIONS	<p>Color & Ratio: Brown 1:1 Total Length: 35.85mm Outlet Type and Size: 2.45mm Material: PP tube + POM element Matched Cartridge: 14ml 1:1,10ml 1:1, 5ml 1:1 syringes Matched oral tip: N3, N4, N8 intra-oral tips</p> <p>Color y proporción: MARRÓN 1: 1 Largo total: 35.85mm Tipo y tamaño de salida: 2.45mm Material: tubo de PP + elemento POM Cartucho combinado: 14ml 1:1,10ml 1:1, 5ml 1:1 jeringa Punta oral combinada: punta intra-oral N3,N4,N8</p>	
PACKAGING	50pcs/bag, 50bags/ctn 50 uds / bolsa, 50 bolsa / ctn	
SHELF LIFE	5 years 5 años	
INTENDED OF USE	For mixing dental cold-whitening cement, dental curing cement, dual-curing core build-up material, dual-curing resin adhesive etc.	

	TECHNICAL DATA SHEET	CODE: TDS_00X
	BROWN MIXING TIPS	DATE: 14/March/20
	PROCLINIC	VERSION: 01

	Para mezclar cemento dental de blanqueamiento en frío, cemento de curado dental, material de acumulación de núcleo de curado dual, adhesivo de resina de curado dual, etc.
SECURITY	Yes Si
TRANSPORT AND STORAGE CONDITIONS	Normal temperature and pressure Temperatura y presión normal