

# Product information



**8690 CSW-FG 660**

<b>Product group</b>	<b>CERASTONE</b>
<b>ISO number</b>	<b>635 314 257 505 025</b>
<b>Single use</b>	<b>no</b>
<b>Medical device class</b>	<b>Ila</b>
<b>Packaging size</b>	<b>10, 100</b>
<b>opt. RPM min</b>	<b>0</b>
<b>opt. RPM max</b>	<b>120000</b>
<b>Weight</b>	<b>0,3 g</b>
<b>Overall length</b>	<b>2 210.000</b>
<b>Custom tariff code</b>	<b>90184910</b>
<b>Country of origin (abbreviation)</b>	<b>CH</b>
<b>Country of Origin (name)</b>	<b>Schweiz</b>
<b>Grit</b>	<b>ultra-fine</b>
<b>Text on label</b>	<b>Arkansas stone for composite processing.</b>
<b>Dimensions (mm)</b>	<b>3 x 8</b>

FB125\_01\_Vorlage\_Produktinformation

EVE Ernst Vetter GmbH  
Neureutstr. 6  
75210 Keltern  
Tel.: +49 7231 97 77 0  
Fax: +49 7231 97 77 99  
Mail: [info@eve-rotary.com](mailto:info@eve-rotary.com)