

The original Tofflemire

Matrix Retainers . Porte-matrices . Matrizenspanner

Portamatrici . Porta-matrices . Porta matriz



Stainless Steel

Junior
+12 Matrix Asstd.

REF

41504

41501

= (41504 + 41528)



Senior
+12 Matrix Asstd.

41505

41502

= (41505 + 41528)




Universal
+12 Matrix Asstd.

41506

41503

= (41506 + 41528)



◊ VASELINE . VASELINA . VASELINA

FOR SMOOTH WORKING . POUR UNE PARFAITE MANIABILITÉ
FÜR EINE ANGENEHME BENÜTZUNG
PER MELIORARE L'IMPIEGO . PARA UNA UTILIZAÇÃO SUELTA
PARA UMA MELHOR UTILIZAÇÃO

Various possibilities of use . Nombreuses possibilités d'emploi
Vermehrte Anwendungsmöglichkeiten . Svariata possibilità di applicazione
Varias modos de empleo . Muitas possibilidades de utilização
Disinfection and sterilization with each use . Désinfection et stérilisation
avant chaque usage . Desinfektion und Sterilisation mit jedem Gebrauch
Desinfección e esterilización con ogni uso . Desinfección y
esterilización . Desinfecção e esterilização com cada uso con cada uso




Tofflemire Matrix Bands

Matrices . Matrizen . Matrici . Matrices . Matrices

Alloy Stainless Steel

Thickness Thin 0.05 mm, x-thin 0.035 mm

Actual sizes



Pack of 12 pieces	REF	REF
Adults	thin	x-thin
# 1	41521	41511
# 2	41522	41512
# 3	41523	41513
# 11	41524	41514
Children		
# 13	41525	41515
# 14	41526	41516
# 15	41527	41517
Assorted		
4 × # 1		
3 × # 2	41528	41518
1 × # 3-15		

 **Produits Dentaires SA**
Rue des Bosquets 18
1800 Vevey . Switzerland



Made in Switzerland
pd-dental.com

For dental professional use . Pour usage dentaire professionnel
Für den professionellen zahnärztlichen Gebrauch . Per uso dentale professionale
Para uso dental profesional . Para uso dental profesional