



Firma / company	Drucker Modell / printer model	V-Print										
		c&b temp	SG clear	splint clear	splint comfort clear	model beige	model 2.0 beige	model fast blue	tray blue	Try-In beige	dentbase pink	cast blue
VOCO	SolFlex 650 / SMP	✓*	✓	✓	✓*	✓	✓*	✓	✓	✓	✓	✓
VOCO	SolFlex 350 / SMP	✓*	✓	✓	✓*	✓	✓*	✓	✓	✓	✓	✓
VOCO	SolFlex 170 / SMP	✓*	✓	✓	✓*	✓	✓*	✓	✓	✓	✓	✓
VOCO	SolFlex 170 HD	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
W2P	SolFlex 650 / PLUS	✓*	✓	✓	✓*	✓	✓*	✓	✓	✓	✓	✓
W2P	SolFlex 350 / PLUS	✓*	✓	✓	✓*	✓	✓*	✓	✓	✓	✓	✓
W2P	SolFlex 170 / PLUS	✓*	✓	✓	✓*	✓	✓*	✓	✓	✓	✓	✓
W2P	SolFlex 163 FullHD	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
ASIGA	MAX UV	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
ASIGA	PRO 4K	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Rapid Shape	D20 II	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Rapid Shape	D30 II	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Rapid Shape	D40 II	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Rapid Shape	D10+	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Rapid Shape	D20+	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Rapid Shape	D30+	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Rapid Shape	D50+	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Ackuretta	SOL	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Microlay	Versus	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Microlay	EVE Pro	—	✓	✓	✓	—	✓	✓	✓	✓	✓	✓
Straumann	P20	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Straumann	P30	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Straumann	P40	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Straumann	P10+	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Straumann	P20+	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Straumann	P30+	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Straumann	P50+	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
MiiCraft	Prime 150Y	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
MiiCraft	Ultra 125Y / 150Y	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
MiiCraft	Hyper 125Y	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
MiiCraft	Profession 250Y	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Shining 3D	AccuFab-L4D	✓	✓	✓	✓	—	✓	—	■	■	✓	✓
Uniz	SLASH 2 PLUS	—	■	✓	✓	—	■	—	—	—	✓	—
Union Tech	EvoDent S300	—	✓	—	✓	—	✓	—	—	—	✓	—

✓ validiert / validation done V-Print Materialsettings sind über Updates seitens des Herstellers verfügbar. / V-Print material settings are available via updates from the manufacturer.

■ in Bearbeitung / in process *Nur mit PowerVat / only with PowerVat

Bitte stellen Sie sicher das Ihre CAM / Slicing-Software in einer aktuellen Version vorliegt. / Please make sure that your CAM / Slicing software is available in a current version.



Firma / company	Drucker Modell / printer model	V-Print										
		c&b temp	SG clear	splint clear	splint comfort clear	model beige	model 2.0 beige	model fast blue	tray blue	Try-In beige	dentbase pink	cast blue
Shera	SHERAprint 20	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Shera	SHERAprint 30	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Shera	SHERAprint 40	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Shera	SHERAprint 20+	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Shera	SHERAprint 30+	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Shera	SHERAprint 50+	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

Dentalwings	D20 II	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Dentalwings	D30 II	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Dentalwings	D40 II	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Dentalwings	D10+	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Dentalwings	D20+	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Dentalwings	D30+	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Dentalwings	D50+	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

DMG	3Delite	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
DMG	3Demax	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

✓ validiert / validation done V-Print Materialsettings müssen vom Anwender beim Druckeranbieter angefragt werden. / V-Print material settings must be requested by the user from the printer supplier.

■ in Bearbeitung / in process *Nur mit PowerVat / only with PowerVat

Bitte stellen Sie sicher das Ihre CAM / Slicing-Software in einer aktuellen Version vorliegt. / Please make sure that your CAM / Slicing software is available in a current version.



TRUSTED PARTNER



Alle Druckpartner finden Sie unter /
All print partners can be found here

www.voco.dental/3dprintingpartners

VOCO GmbH
Anton-Flettner-Str. 1-3
27472 Cuxhaven
Germany

Phone +49 (4721) 719-0
Fax +49 (4721) 719-140
e-mail: marketing@voco.com
www.voco.dental

