



## FICHA TÉCNICA

<b>CERA ESQUELÉTICOS REJILLA PROCLINIC</b>		
<b>REFERENCIAS</b>	<b>H20072:</b> CUADRADA <b>H20074:</b> REDONDA <b>H20076:</b> REDONDA 40 U	
<b>FABRICANTE / RESPONSABLE</b>	<b>Al Dente Dentalprodukte GmbH</b> Borsigstr. 1 - D-38644 Goslar +49 (0) 5321-80031 info@aldente.de Emergency Telephone Number: +49 (0) 171 8637790	
<b>CLASIFICACIÓN</b>	Producto clase protésico.	
<b>DESCRIPCIÓN</b>	Mallas de retención de cera.	
<b>COMPOSICIÓN</b>	Parafina.	
<b>PRESENTACIÓN</b>	H20072-4: 20 u por envase. H20076: 40 u por envase.	
<b>FINALIDAD</b>	Rejillas de cera para retención en técnica de esqueléticos.	
<b>INSTRUCCIONES DE USO</b>	Ver instrucciones en el producto.	
<b>ALMACENAMIENTO Y TRANSPORTE</b>	No se requieren condiciones especiales de almacenamiento y transporte.	

## TECHNICAL DATA SHEET

<b>WAX RETENTIONS MESHES PROCLINIC</b>		
<b>REFERENCES</b>	<b>H20072: SQUARE</b> <b>H20074: ROUND</b> <b>H20076: ROUND 40 U</b>	
<b>MANUFACTURER</b>	<b>Al Dente Dentalprodukte GmbH</b> Borsigstr. 1 - D-38644 Goslar +49 (0) 5321-80031 info@aldente.de Emergency Telephone Number: +49 (0) 171 8637790	
<b>CLASSIFICATION</b>	Prosthetic class product.	
<b>DESCRIPTION</b>	Wax retention meshes.	
<b>COMPOSITION</b>	Paraffin.	
<b>PACKAGING</b>	H20072-4: 20 u per container. H20076: 40 u per container.	
<b>INTENDED USE</b>	Wax grids for retention in the skeletal technique.	
<b>INSTRUCTIONS FOR USE</b>	See instructions on the product.	
<b>HANDLING AND STORAGE CONDITIONS</b>	No special storage and transport conditions are required.	