



## FICHA TÉCNICA

<b>CERA ESQUELETICOS GANCHOS PROCLINIC</b>		
<b>REFERENCIAS</b>	<b>H20078:</b> MOLARES <b>H20080:</b> PREMOLARES <b>H20082:</b> BONYHARD <b>H20084:</b> ANULARES RECTO	
<b>FABRICANTE / RESPONSABLE</b>	<b>Al Dente Dentalprodukte GmbH</b> Borsigstr. 1 - D-38644 Goslar +49 (0) 5321-80031 info@aldente.de Emergency Telephone Number: +49 (0) 171 8637790	
<b>CLASIFICACIÓN</b>	Producto clase protésico.	
<b>DESCRIPCIÓN</b>	Cera calcinable adaptable al modelo sin tensiones. La transparencia de la cera permite ver marcas a través de ella.	
<b>COMPOSICIÓN</b>	Parafina.	
<b>PRESENTACIÓN</b>	200 u por envase.	
<b>FINALIDAD</b>	Agentes auxiliares para la fabricación en la tecnología dental.	
<b>INSTRUCCIONES DE USO</b>	Ver instrucciones en el producto.	
<b>ALMACENAMIENTO Y TRANSPORTE</b>	No se requieren condiciones especiales de almacenamiento y transporte.	

## TECHNICAL DATA SHEET

<b>SKELETON WAX HOOKS PROCLINIC</b>		
<b>REFERENCES</b>	<b>H20078:</b> MOLARS <b>H20080:</b> PREMOLARS <b>H20082:</b> BONYHARD <b>H20084:</b> RING STRAIGHT	
<b>MANUFACTURER</b>	<b>Al Dente Dentalprodukte GmbH</b> Borsigstr. 1 - D-38644 Goslar +49 (0) 5321-80031 info@aldente.de Emergency Telephone Number: +49 (0) 171 8637790	
<b>CLASSIFICATION</b>	Prosthetic class product.	
<b>DESCRIPTION</b>	Castable wax adaptable to the model without tensions. The transparency of the wax allows you to see marks through it.	
<b>COMPOSITION</b>	Paraffin.	
<b>PACKAGING</b>	200 u per-container.	
<b>INTENDED USE</b>	Auxiliary agents for manufacturing in dental technology.	
<b>INSTRUCTIONS FOR USE</b>	See instructions on the product.	
<b>HANDLING AND STORAGE CONDITIONS</b>	No special storage and transport conditions are required.	