

TruCast[™]

Low Profile Buccal Tubes

- Single Non-Convertible Bondable
- •Lower profile design to avoid occlusal interference
- Chamfered slot for easy archwire placement
- •Streamlined design for maximum patient comfort
- Color coded scribe line for height and angulation reference point.
 Red=upper, Yellow=lower
- •Bondable are indented 80 gauge mesh pad for an anatomically formed fit and maximum bond strength
- Occlusal and gingival positioning notches for ease of placement and handling
- •Comfort designed ball hook enhances patient comfort and avoids gingival impingement
- Mesial/distal dimension 3.6 mm

Positioning Notch To assist with placement and avoid interference

Removable Color Coded Identification For quick and accurate placemen



Low Profile Hook For patient comfort Indented Mesh Pad For a precise contoured fit

Lower Profile Design To reduce occlusal interference

Mesial/Distal Dimension 3.6 mm

Chamfered Entry Guides wires into buccal tubes