

## FICHA TÉCNICA

|                             |                                                                             |                                                                                    |
|-----------------------------|-----------------------------------------------------------------------------|------------------------------------------------------------------------------------|
| TAZA DE ALGINATO PROCLINIC  |                                                                             |  |
| REFERENCIA                  | 0111                                                                        |                                                                                    |
| FABRICANTE / RESPONSABLE    | LARIDENT SRL<br>Via Lamaneigra 12R - 16030 Tribogna (GE)<br>+39 0185 939481 |                                                                                    |
| CLASIFICACIÓN               | Producto no Sanitario (Clase consumo).                                      |                                                                                    |
| DESCRIPCIÓN                 | Taza flexible para alginatos.                                               |                                                                                    |
| ESPECIFICACIONES            | Color: verde<br>Sin látex<br>Diámetro: 12,5 cm<br>Altura: 9 cm              |                                                                                    |
| COMPOSICIÓN                 | PVC PT65 No tóxico.                                                         |                                                                                    |
| EMBALAJE                    | 1 unidad.                                                                   |                                                                                    |
| PROPIEDADES TÉRMICAS        | A 30° C tiende a ablandarse.                                                |                                                                                    |
| PROPIEDADES QUÍMICAS        | Resiste todos los productos químicos de limpieza estándar.                  |                                                                                    |
| ALMACENAMIENTO Y TRANSPORTE | No requiere condiciones especiales de almacenamiento y transporte.          |                                                                                    |

## TECHNICAL DATA SHEET

|                         |                                                                             |                                                                                    |
|-------------------------|-----------------------------------------------------------------------------|------------------------------------------------------------------------------------|
| ALGINATE BOWL PROCLINIC |                                                                             |  |
| REFERENCE               | 0111                                                                        |                                                                                    |
| MANUFACTURER            | LARIDENT SRL<br>Via Lamaneigra 12R - 16030 Tribogna (GE)<br>+39 0185 939481 |                                                                                    |
| CLASSIFICATION          | Non-Sanitary Product (Consumer Class).                                      |                                                                                    |
| DESCRIPTION             | Flexible bowl for alginates.                                                |                                                                                    |
| SPECIFICATIONS          | Color: green<br>Latex Free<br>Diameter: 12,5 cm<br>Height: 9 cm             |                                                                                    |
| COMPOSITION             | PVC PT65 Non-toxic                                                          |                                                                                    |
| PACKAGING               | 1 unit.                                                                     |                                                                                    |
| THERMAL PROPERTIES      | At 30° C it tends to soften.                                                |                                                                                    |
| CHEMICAL PROPERTIES     | It resists all standard chemical products for cleaning.                     |                                                                                    |
| MA                      | It does not require special storage and transportation conditions.          |                                                                                    |