

FICHA TÉCNICA

YESO DE ORTODONCIA EXTRABLANCO TIPO III/3		
REFERENCIAS	01148	
FABRICANTE / RESPONSABLE	MADESPA, S.A. C/ Río Jarama, 120, 45007, Toledo España	
CLASIFICACIÓN	Producto No Sanitario (Clase protésico).	
DESCRIPCIÓN	Yeso producido con roca seleccionada con el fin de garantizar un producto blanquísimo, homogéneo y compacto. El modelo no se deforma durante la polimerización de las resinas acrílicas. Líder en el mercado por su fina granulometría y su color blanco brillante. Conforme norma I.S.O. 6873.	
ESPECIFICACIÓN	Proporción de mezcla agua/polvo: 26 ml/100 gr Tiempo de trabajo: 6'-8' Expansión de fraguado (2h): ≤0,20 % Resistencia a la compresión (1h): ~33 MPa	
COMPOSICIÓN	Sulfato de calcio.	
PRESENTACIÓN	1 unidad de 5 kg.	
VIDA ÚTIL	3 años.	
INSTRUCCIONES DE USO	Poner agua en una taza y añadir lentamente el polvo según la proporción recomendada. Se recomienda utilizar vibrador.	
FINALIDAD PREVISTA	Yeso duro de tipo 3 para la realización de modelos de ortodoncia, color blanco brillante.	
ALMACENAMIENTO Y TRANSPORTE	Conservar los envases en un lugar seco, bien ventilado y protegido de la humedad.	

TECHNICAL DATA SHEET

EXTRA-WHITE HARD GYPSUM TYPE III/3		
REFERENCES	01148	
MANUFACTURER	MADESPA, S.A. C/ Río Jarama, 120, 45007, Toledo Spain	
CLASSIFICATION	Non-Sanitary Product (Prosthetic Class).	
DESCRIPTION	Gypsum produced with selected rock in order to guarantee a very white, homogeneous and compact product. The model does not deform during the polymerization of acrylic resins. Market leader for its fine granulometry and bright white color. Conforms to I.S.O. 6873 standard.	
SPECIFICATION	Mixing ratio water/powder: 26 ml/100 gr. Working time: 6'-8'. Setting expansion (2h): ≤0,20 %. Compressive strength (1h): ~33 MPa	
COMPOSITION	Calcium sulfate.	
PACKAGE	1 unit of 5kg.	
SHELF LIFE	3 years.	
INSTRUCTIONS FOR USE	Put water in a cup and slowly add the powder according to the recommended proportion. The use of a vibrator is recommended.	
INTENDED USE	Type 3 hard plaster for the realization of orthodontic models, bright white color.	
STORAGE AND TRANSPORT	Store containers in a dry, well-ventilated place, protected from humidity.	