

FICHA TÉCNICA

DISPENSADOR PELLETS DE ALGODÓN PROCLINIC		
REFERENCIA	0115	
FABRICANTE / RESPONSABLE	LARIDENT SRL Via Lamaneigra 12R - 16030 Tribogna (GE) +39 0185 939481 www.larident.it	
CLASIFICACIÓN	Producto no Sanitario (Clase consumo).	
DESCRIPCIÓN	Dispensador pellets de algodón.	
ESPECIFICACIONES	Medidas: 5x5 cm h 5.3. Color Blanco	
COMPOSICIÓN	Fabricado en plástico autoclavable. Marco: POLIAMIDA. Base: CAUCHO TERMOPLÁSTICO. Muelle y placa de acero en poliamida.	
PRESENTACIÓN	1 unidad.	
PROPIEDADES TÉRMICAS	Autoclavable a 121°C.	
PROPIEDADES QUÍMICAS	Resiste a todos los esterilizadores químicos.	
ALMACENAMIENTO Y TRANSPORTE	No requiere condiciones especiales de almacenamiento y transporte.	

TECHNICAL DATA SHEET

PROCLINIC COTTON PELLETS DISPENSER		
REFERENCE	0115	
MANUFACTURER	LARIDENT SRL Via Lamaneigra 12R - 16030 Tribogna (GE) +39 0185 939481 www.larident.it	
CLASSIFICATION	Non-Sanitary Product (Consumer Class).	
DESCRIPTION	Dispenser pellets Cotton.	
SPECIFICATIONS	Measurements: cm 5x5 h 5.3 Color: White	
COMPOSITION	Frame: polyamide. Base: thermoplastic rubber. Steel spring and plate in polyamide.	
PACKAGING	1 unit.	
THERMAL PROPERTIES	Autoclavable at 121°C.	
CHEMICAL PROPERTIES	It resists all chemical sterilizers.	
STORAGE AND TRANSPORT	It does not require special storage and transportation conditions.	