

FICHA TÉCNICA

CEPILLOS DE NYLON PROCLINIC		
REFERENCIAS	01201: Duros Azules. 01202: Medios Púrpura. 01203: Suaves Rosas.	
FABRICANTE / RESPONSABLE	STODDARD MANUFACTURING Co. Ltd. Blackhorse road, letchworth, herts, sg6 1hb, England - Tel: + 44 (0) 1462 686221 www.stoddard.co.uk	
CLASIFICACIÓN	Producto Sanitario Clase IIa, Directiva 93/42/CEE.	
DESCRIPCIÓN	- Cepillos para contra-ángulo fabricados en nylon de alta calidad con diferentes grados de dureza para adaptarse a todas las necesidades clínicas. - Color azul para una mayor dureza, púrpura para una dureza intermedia y rosa para una dureza suave. - Perfectos para el uso con pasta de pulido. - Libres de látex.	
ESPECIFICACIÓN	<u>Color:</u> Azul, Púrpura, Rosas.	
COMPOSICIÓN	Filamento de nailon - Polihexametileno dodecanamida, vástago de latón, virola de latón, anillo de latón.	
PRESENTACIÓN	100 unidades por envase.	
FINALIDAD PREVISTA	Se utilizan para la limpieza profiláctica o el pulido de dientes naturales y superficies de obturación en la boca. Autoclavar mediante un procedimiento validado a 134°C.	
ALMACENAMIENTO Y TRANSPORTE	No se requieren condiciones especiales de almacenamiento y transporte.	

TECHNICAL DATA SHEET

NYLON BRUSHES PROCLINIC		
REFERENCES	01201: Hard Blue. 01212: Medium Purple. 01203: Soft roses.	
MANUFACTURER	STODDARD MANUFACTURING Co. Ltd. Blackhorse road, letchworth, herts, sg6 1hb, England - Tel: + 44 (0) 1462 686221 www.stoddard.co.uk	
CLASSIFICATION	Medical Device, Class IIa, Directive 93/42/EEC.	
DESCRIPTION	- Contra-angle brushes made of high quality nylon with different degrees of hardness to adapt to all clinical needs. - Blue for a higher hardness, purple for an intermediate hardness and pink for a soft hardness. - Perfect for use with polishing paste. - Latex free.	
SPECIFICATION	<u>Color:</u> Blue, Purple, Roses.	
COMPOSITION	Nylon filament: Polyhexamethylene dodecanamide, brass stem, brass ferrule, brass ring.	
PACKAGE	100 units per-container.	
INTENDED USE	Used for the prophy cleaning or polishing of natural teeth and filled surfaces in the mouth. Autoclave using a validated procedure at 134°C.	
STORAGE AND TRANSPORT	No special storage and transport conditions are required.	