

FICHA TÉCNICA

COPAS DE GOMA PROCLINIC		
REFERENCIAS	01261: Suaves rosas. 01262: Medias púpura. 01263: Duras azules.	
FABRICANTE / RESPONSABLE	STODDARD BLACKHORSE ROAD, LETCHWORTH, HERTS, SG6 1HB, ENGLAND Tel: + 44 (0) 1462 686221 www.stoddard.co.uk	
CLASIFICACIÓN	Producto Sanitario Clase IIa, Directiva 93/42/CEE.	
DESCRIPCIÓN	- Copas de goma para contra-ángulo fabricadas en 3 tipos diferentes de polímero sintético que le aportan distintos grados de dureza para adaptarse a cada necesidad clínica. - Color azul para una mayor dureza, púpura para una dureza intermedia y rosa para una dureza blanda. - Perfectas para el uso con pasta de pulido. - Libres de látex. - Velocidad máxima 10.000 rpm	
COMPOSICIÓN	Elastómero termoplástico, vástago de latón, polipropileno, anillo de caucho de silicona.	
PRESENTACIÓN	36 unidades por envase.	
FINALIDAD PREVISTA	Son productos de un solo uso para el pulido profiláctico de dientes naturales, oro y empastes de amalgama. Autoclavar mediante un procedimiento validado a 134 ⁰ C.	
ALMACENAMIENTO Y TRANSPORTE	No se requieren condiciones especiales de almacenamiento y transporte.	

TECHNICAL DATA SHEET

RUBBER CUPS PROCLINIC			
REFERENCES	01261: Soft roses. 01262: Medium Purple. 01263: Hard blue.		
MANUFACTURER	STODDARD BLACKHORSE ROAD, LETCHWORTH, HERTS, SG6 1HB, ENGLAND Tel: + 44 (0) 1462 686221 www.stoddard.co.uk		
CLASSIFICATION	Class IIa Medical Device, Directive 93/42/EEC.		
DESCRIPTION	-Rubber cups for counter-manufactured in 3 different types of synthetic polymer that provide different degrees of hardness to adapt to each clinical need. - Blue for a higher hardness, purple for an intermediate hardness and pink for a soft hardness. - Perfect for use with polishing paste. - Latex free. - Maximum speed 10,000 rpm		
COMPOSITION	Thermoplastic elastomer, brass stem, polypropylene, silicone rubber ring.		
PACKAGE	36 units per-container.		
INTENDED USE	They are single-use products for the prophylactic polishing of natural teeth, gold and amalgam fillings. Autoclave by a validated procedure at 134°C.		
HANDLING AND STORAGE CONDITIONS	No special storage and transport conditions are required.		