

FICHA TÉCNICA

PULIDORES PROCLINIC		
REFERENCIAS	0180: DESBASTAR. TAZA NEGRA. 0181: DESBASTAR. PUNTA NEGRA. 0184: PULIDO. MINIPUNTA MARRÓ. 0182: PULIDO. TAZA MARRÓN. 0183: PULIDO. PUNTA MARRÓN. 0187: BRILLO. FINAL MINIPUNTA VERDE. 0185: BRILLO. FINAL TAZA VERDE. 0186: BRILLO. FINAL PUNTA VERDE.	
FABRICANTE / RESPONSABLE	EVE Ernst Vetter GmbH Neureutstr. 6- 75210 Keltern Germany Tel.: +49 7231 97 77 -0 Fax: +49 7231 97 77 -99 E-Mail: info@eve-rotary.com	
CLASIFICACIÓN	Producto Sanitario Clase IIa, bajo reglamento (EU) 2017/745.	
DESCRIPCIÓN	Pulidores montados en silicona, para el acabado de amalgamas, oro, etc.	
ESPECIFICACIÓN	Tamaños: 0180: DESBASTAR. TAZA NEGRA: 6.5 x 9.5 mm 0181: DESBASTAR. PUNTA NEGRA: 5 x 10 mm 0184: PULIDO. MINIPUNTA MARRÓ: 3 x 7.5 mm 0182: PULIDO. TAZA MARRÓN: 6.5 x 9.5 mm 0183: PULIDO. PUNTA MARRÓN: 5 x 10 mm 0187: BRILLO. FINAL MINIPUNTA VERDE: 3 x 7.5 mm 0185: BRILLO. FINAL TAZA VERDE: 6.5 x 9.5 mm 0186: BRILLO. FINAL PUNTA VERDE: 5 x 10 mm	
PRESENTACIÓN	10 unidades por caja	
FINALIDAD	Para pulir oro, amalgamas y aleaciones no preciosas.	
ALMACENAMIENTO Y TRANSPORTE	No requiere condiciones especiales de almacenamiento y transporte.	

TECHNICAL DATA SHEET

POLISHERS PROCLINIC			
REFERENCES	0180: TOUGH. BLACK MUG 0181: TOUGH. BLACK POINT 0184: POLISHED. BROWN MINI POINT 0182: POLISHED. BROWN MUG 0183: POLISHED. BROWN POINT 0187: BRIGHTNESS. FINAL GREEN MINI POINT 0185: BRIGHTNESS. FINAL GREEN CUP 0186: BRIGHTNESS. FINAL GREEN POINT		
MANUFACTURER	EVE Ernst Vetter GmbH Neureutstr. 6- 75210 Keltern Germany Tel.: +49 7231 97 77 -0 Fax: +49 7231 97 77 -99 E-Mail: info@eve-rotary.com		
CLASSIFICATION	Class IIa Medical Device, under regulation (EU) 2017/745.		
DESCRIPTION	Silicone mounted polishers, for finishing amalgams, gold, etc.		
SPECIFICATIONS	Sizes: 0180: TOUGH. BLACK CUP: 6.5 x 9.5mm 0181: TOUGH. BLACK TIP: 5 x 10 mm 0184: POLISHED. BROWN MINI NIP: 3 x 7.5 mm 0182: POLISHED. BROWN CUP: 6.5 x 9.5mm 0183: POLISHED. BROWN TIP: 5 x 10 mm 0187: BRIGHTNESS. FINAL GREEN MINI NIP: 3 x 7.5 mm 0185: BRIGHTNESS. FINAL GREEN CUP: 6.5 x 9.5mm 0186: BRIGHTNESS. GREEN TIP FINAL: 5 x 10 mm		
PACKAGING	10 units per package.		
INTENDED USE	For polishing gold, amalgam and non-precious alloys.		
HANDLING AND STORAGE CONDITIONS	It does not require special storage and transport conditions.		