

FICHA TÉCNICA

PAPEL DE ARTICULAR FINO PROCLINIC		
REFERENCIAS	01954: rojo, grosor 40 µm. 01955: azul, grosor 40 µm.	
FABRICANTE / RESPONSABLE	Coltène/Whaledent GmbH + Co. KG Raiffeisenstraße 30 89129 Langenau/Alemania	
CLASIFICACIÓN	Producto Sanitario Clase I, según Reglamento 2017/745/UE.	
DESCRIPCIÓN DEL PRODUCTO	El papel de articular ultrafino para el control y verificación final de la oclusión marca Proclinic, permite visualizar de forma precisa los contactos prematuros.	
ESPECIFICACIONES	- <u>Dimensiones:</u> 40 µ de espesor 22 mm x 110 mm o 18 mm x 50 mm - <u>Aspecto:</u> Bordes de corte lisos Apenas sobresalen fibras Intensidad de color homogéneo Espiral homogénea - <u>Color:</u> rojo o azul.	
ENVASE	Envases conteniendo 120 unidades (40 micras).	
VIDA ÚTIL	3 años	
FINALIDAD PREVISTA	Marcar los contactos dentales en las restauraciones dentales para visualizar interferencias en las fuerzas masticatorias.	
INFORMACIÓN SOBRE LA SEGURIDAD	Producto de uso por los profesionales dentales. Este producto es bien conocido por el usuario y es seguro.	
ALMACENAMIENTO Y TRANSPORTE	Conservar por debajo de 25°C. Proteger de la luz solar directa y de la humedad.	

TECHNICAL DATA SHEET

PROCLINIC THIN ARTICULATING PAPER		
REFERENCES	01954: red, thickness 40 µm. 01955: blue, thickness 40 µm.	
MANUFACTURER	Coltène/Whaledent GmbH + Co. KG Raiffeisenstraße 30 89129 Langenau/Alemania	
CLASSIFICATION	Class I Medical Device, according to Regulation 2017/745/EU.	
DESCRIPTION	The Proclinic ultra-fine articulating paper for the final check and verification of the occlusion allows precise visualisation of premature contacts.	
SPECIFICATION	- Dimensions: 40 µ thickness 22 mm x 110 mm or 18 mm x 50 mm - Appearance: - Smooth cut edges - Hardly protruding fibers - Homogeneous color intensity - Homogeneous spiral - Color: red or blue.	
PRESENTATION	Packs containing 120 units (40 microns).	
SHELF LIFE	3 years.	
INTENDED USE	Mark tooth contacts on dental restorations to visualise interferences in masticatory forces.	
SECURITY INFORMATION	Product for use by dental professionals. This product is well known to the user and is safe.	
STORAGE AND TRANSPORT	Store below 25°C. Protect from direct sunlight and moisture.	