

FICHA TÉCNICA

PUNTAS DE GUTTAPERCHA BESTDENT																																																																																			
REFERENCIAS	Ver anexo.																																																																																		
FABRICANTE / RESPONSABLE	Shanghai Xingyu Medical Equipment Co., Ltd. 3 Floor, Building 21, No. 3825, Xin Zhuan Road 201601 Shanghai, China																																																																																		
CLASIFICACIÓN	Producto Sanitario clase IIa, conforme la Directiva 93/42/CEE.																																																																																		
DESCRIPCIÓN	Puntas de guttapercha Iso, utilizadas para la obturación del canal radicular en el tratamiento de endodoncia. Código de colores ISO para una identificación fácil.																																																																																		
ESPECIFICACIONES	<table><tr><th>Tamaño</th><th>Diámetro</th><th>Color</th><th>Tamaño</th><th>Diámetro</th><th>Color</th></tr><tr><td>15</td><td>0.15</td><td>Blanco</td><td>45</td><td>0.45</td><td>Blanco</td></tr><tr><td>20</td><td>0.20</td><td>Amarillo</td><td>50</td><td>0.50</td><td>Amarillo</td></tr><tr><td>25</td><td>0.25</td><td>Rojo</td><td>55</td><td>0.55</td><td>Rojo</td></tr><tr><td>30</td><td>0.30</td><td>Azul</td><td>60</td><td>0.60</td><td>Azul</td></tr><tr><td>35</td><td>0.35</td><td>Verde</td><td>70</td><td>0.70</td><td>Verde</td></tr><tr><td>40</td><td>0.40</td><td>Negro</td><td>80</td><td>0.80</td><td>Negro</td></tr></table> Longitud: 28±2mm Diámetro significa el diámetro de la punta del producto. (Taper 0.02) <table><tr><th>Standard specification</th><th>15</th><th>20</th><th>25</th><th>30</th><th>35</th><th>40</th><th>45</th><th>50</th><th>55</th><th>60</th><th>70</th><th>80</th></tr><tr><td>Length mm</td><td>28</td><td>28</td><td>28</td><td>28</td><td>28</td><td>28</td><td>28</td><td>28</td><td>28</td><td>28</td><td>28</td><td>28</td></tr><tr><td>Tip diameter mm</td><td>0.15</td><td>0.20</td><td>0.25</td><td>0.30</td><td>0.35</td><td>0.40</td><td>0.45</td><td>0.50</td><td>0.55</td><td>0.60</td><td>0.70</td><td>0.80</td></tr></table>		Tamaño	Diámetro	Color	Tamaño	Diámetro	Color	15	0.15	Blanco	45	0.45	Blanco	20	0.20	Amarillo	50	0.50	Amarillo	25	0.25	Rojo	55	0.55	Rojo	30	0.30	Azul	60	0.60	Azul	35	0.35	Verde	70	0.70	Verde	40	0.40	Negro	80	0.80	Negro	Standard specification	15	20	25	30	35	40	45	50	55	60	70	80	Length mm	28	28	28	28	28	28	28	28	28	28	28	28	Tip diameter mm	0.15	0.20	0.25	0.30	0.35	0.40	0.45	0.50	0.55	0.60	0.70	0.80
Tamaño	Diámetro	Color	Tamaño	Diámetro	Color																																																																														
15	0.15	Blanco	45	0.45	Blanco																																																																														
20	0.20	Amarillo	50	0.50	Amarillo																																																																														
25	0.25	Rojo	55	0.55	Rojo																																																																														
30	0.30	Azul	60	0.60	Azul																																																																														
35	0.35	Verde	70	0.70	Verde																																																																														
40	0.40	Negro	80	0.80	Negro																																																																														
Standard specification	15	20	25	30	35	40	45	50	55	60	70	80																																																																							
Length mm	28	28	28	28	28	28	28	28	28	28	28	28																																																																							
Tip diameter mm	0.15	0.20	0.25	0.30	0.35	0.40	0.45	0.50	0.55	0.60	0.70	0.80																																																																							
COMPOSICIÓN	Resina natural, Óxido de Zinc, Sulfato de Bario, Colorante.																																																																																		
ENVASE	Unidades Variadas.																																																																																		
VIDA ÚTIL	3 años.																																																																																		
FINALIDAD PREVISTA	Sólo para uso dental y para la obturación del canal radicular en el tratamiento de endodoncia.																																																																																		
ALMACENAMIENTO Y TRANSPORTE	No requiere condiciones especiales de almacenamiento y transporte.																																																																																		

TECHNICAL DATA SHEET

BESTDENT GUTTA-PERCHA TIPS																																																																																			
REFERENCES	See annex 1.																																																																																		
MANUFACTURER	Shanghai Xingyu Medical Equipment Co., Ltd. 3 Floor, Building 21, No. 3825, Xin Zhuan Road 201601 Shanghai, China																																																																																		
CLASSIFICATION	Medical device class IIa, according to Directive 93/42/EEC.																																																																																		
DESCRIPTION	Iso gutta-percha tips, used for root canal obturation in endodontic treatment. ISO colour coded for easy identification.																																																																																		
SPECIFICATION	<table><thead><tr><th>Tamaño</th><th>Diámetro</th><th>Color</th><th>Tamaño</th><th>Diámetro</th><th>Color</th></tr></thead><tbody><tr><td>15</td><td>0.15</td><td>Blanco</td><td>45</td><td>0.45</td><td>Blanco</td></tr><tr><td>20</td><td>0.20</td><td>Amarillo</td><td>50</td><td>0.50</td><td>Amarillo</td></tr><tr><td>25</td><td>0.25</td><td>Rojo</td><td>55</td><td>0.55</td><td>Rojo</td></tr><tr><td>30</td><td>0.30</td><td>Azul</td><td>60</td><td>0.60</td><td>Azul</td></tr><tr><td>35</td><td>0.35</td><td>Verde</td><td>70</td><td>0.70</td><td>Verde</td></tr><tr><td>40</td><td>0.40</td><td>Negro</td><td>80</td><td>0.80</td><td>Negro</td></tr></tbody></table> Length: 28-+2mm. Diameter means the diameter of the tip of the product. (Taper 0.02) <table><thead><tr><th>Standard specification</th><th>15</th><th>20</th><th>25</th><th>30</th><th>35</th><th>40</th><th>45</th><th>50</th><th>55</th><th>60</th><th>70</th><th>80</th></tr></thead><tbody><tr><td>Length mm</td><td>28</td><td>28</td><td>28</td><td>28</td><td>28</td><td>28</td><td>28</td><td>28</td><td>28</td><td>28</td><td>28</td><td>28</td></tr><tr><td>Tip diameter mm</td><td>0.15</td><td>0.20</td><td>0.25</td><td>0.30</td><td>0.35</td><td>0.40</td><td>0.45</td><td>0.50</td><td>0.55</td><td>0.60</td><td>0.70</td><td>0.80</td></tr></tbody></table>		Tamaño	Diámetro	Color	Tamaño	Diámetro	Color	15	0.15	Blanco	45	0.45	Blanco	20	0.20	Amarillo	50	0.50	Amarillo	25	0.25	Rojo	55	0.55	Rojo	30	0.30	Azul	60	0.60	Azul	35	0.35	Verde	70	0.70	Verde	40	0.40	Negro	80	0.80	Negro	Standard specification	15	20	25	30	35	40	45	50	55	60	70	80	Length mm	28	28	28	28	28	28	28	28	28	28	28	28	Tip diameter mm	0.15	0.20	0.25	0.30	0.35	0.40	0.45	0.50	0.55	0.60	0.70	0.80
Tamaño	Diámetro	Color	Tamaño	Diámetro	Color																																																																														
15	0.15	Blanco	45	0.45	Blanco																																																																														
20	0.20	Amarillo	50	0.50	Amarillo																																																																														
25	0.25	Rojo	55	0.55	Rojo																																																																														
30	0.30	Azul	60	0.60	Azul																																																																														
35	0.35	Verde	70	0.70	Verde																																																																														
40	0.40	Negro	80	0.80	Negro																																																																														
Standard specification	15	20	25	30	35	40	45	50	55	60	70	80																																																																							
Length mm	28	28	28	28	28	28	28	28	28	28	28	28																																																																							
Tip diameter mm	0.15	0.20	0.25	0.30	0.35	0.40	0.45	0.50	0.55	0.60	0.70	0.80																																																																							
COMPOSITION	Natural resin, Zinc oxide, Barium sulphate, Colouring agent.																																																																																		
PACKAGING	Various units.																																																																																		
SHELF LIFE	3 years.																																																																																		
INTENDED USE	For dental use only and for root canal filling in endodontic treatment.																																																																																		
STORAGE AND TRANSPORT	No special storage and transport conditions required.																																																																																		

ANEXO/ ANNEX 1:

Referencia	Descripción
100400	GUTTAPERCHA CONICIDAD 04 N.15 BESTDENT
100401	GUTTAPERCHA CONICIDAD 04 N.25 BESTDENT
100402	GUTTAPERCHA CONICIDAD 04 N.30 BESTDENT
100403	GUTTAPERCHA CONICIDAD 04 N.35 BESTDENT
100404	GUTTAPERCHA CONICIDAD 04 N.40 BESTDENT
100405	GUTTAPERCHA CONIC. 04 N.15 a 40 BESTDENT
100406	GUTTAPERCHA CONICIDAD 04 N.45 BESTDENT
100407	GUTTAPERCHA CONIC. 04 N.45 a 80 BESTDENT
100408	GUTTAPERCHA CONICIDAD 06 N.15 BESTDENT
100409	GUTTAPERCHA CONICIDAD 06 N.20 BESTDENT
100410	GUTTAPERCHA CONICIDAD 06 N.25 BESTDENT
100411	GUTTAPERCHA CONICIDAD 06 N.30 BESTDENT
100412	GUTTAPERCHA CONICIDAD 06 N.35 BESTDENT
100413	GUTTAPERCHA CONICIDAD 06 N.40 BESTDENT
100414	GUTTAPERCHA CONIC. 06 N.15 a 40 BESTDENT
100415	GUTTAPERCHA TAPER F1 BESTDENT
100416	GUTTAPERCHA TAPER F2 BESTDENT
100417	GUTTAPERCHA TAPER F3 BESTDENT
100418	GUTTAPERCHA TAPER SURT F1-F2-F3 BD
100419	GUTTAPERCHA TAPER SURT F4-F5 BESTDENT
100420	GUTTAPERCHA N.15 BESTDENT
100421	GUTTAPERCHA N.20 BESTDENT
100422	GUTTAPERCHA N.25 BESTDENT
100423	GUTTAPERCHA N.30 BESTDENT
100424	GUTTAPERCHA N.35 BESTDENT
100425	GUTTAPERCHA N.40 BESTDENT
100426	GUTTAPERCHA SURTIDA 15-40 BESTDENT
100427	GUTTAPERCHA N.45 BESTDENT
100428	GUTTAPERCHA N.50 BESTDENT
100429	GUTTAPERCHA N.55 BESTDENT
100430	GUTTAPERCHA N.60 BESTDENT
100431	GUTTAPERCHA N.70 BESTDENT
100432	GUTTAPERCHA N.80 BESTDENT
100433	GUTTAPERCHA SURTIDA 45-80 BESTDENT
100434	GUTTAPERCHA BESTDENT -XF-
100435	GUTTAPERCHA BESTDENT -F-
100436	GUTTAPERCHA BESTDENT -M-
100464	GUTTAPERCHA REC PARA R25 (60u) BD
100465	GUTTAPERCHA REC PARA R40 (60u) BD
100466	GUTTAPERCHA REC PARA R50 (60u) BD
100467	GUTTAPERCHA REC SURTIDA (60u) BD
100472	GUTTAPERCHA CONICIDAD 04 N.20 BESTDENT