

FICHA TÉCNICA

PUNTAS DE PAPEL BESTDENT																																																																																			
REFERENCIAS	Ver anexo 1.																																																																																		
FABRICANTE / RESPONSABLE	Shanghai Xingyu Medical Equipment Co., Ltd. 3 Floor, Building 21, No. 3825, Xin Zhuan Road 201601 Shanghai, China																																																																																		
CLASIFICACIÓN	Producto Sanitario clase IIa, conforme la Directiva 93/42/CEE.																																																																																		
DESCRIPCIÓN	Puntas de papel absorbentes para secar los canales radiculares preparados con instrumentos convencionales antes de su obturación. Código de colores ISO para una identificación fácil.																																																																																		
ESPECIFICACIONES	<table><tr><td>Tamaño</td><td>Diámetro</td><td>Color</td><td>Tamaño</td><td>Diámetro</td><td>Color</td></tr><tr><td>15</td><td>0.15</td><td>Blanco</td><td>45</td><td>0.45</td><td>Blanco</td></tr><tr><td>20</td><td>0.20</td><td>Amarillo</td><td>50</td><td>0.50</td><td>Amarillo</td></tr><tr><td>25</td><td>0.25</td><td>Rojo</td><td>55</td><td>0.55</td><td>Rojo</td></tr><tr><td>30</td><td>0.30</td><td>Azul</td><td>60</td><td>0.60</td><td>Azul</td></tr><tr><td>35</td><td>0.35</td><td>Verde</td><td>70</td><td>0.70</td><td>Verde</td></tr><tr><td>40</td><td>0.40</td><td>Negro</td><td>80</td><td>0.80</td><td>Negro</td></tr></table> Longitud: 28±2mm Diámetro significa el diámetro de la punta del producto. (Taper 0.02, Taper 0.04, Taper 0.6) <table><tr><td>Standard specification</td><td>15</td><td>20</td><td>25</td><td>30</td><td>35</td><td>40</td><td>45</td><td>50</td><td>55</td><td>60</td><td>70</td><td>80</td></tr><tr><td>Length mm</td><td>28</td><td>28</td><td>28</td><td>28</td><td>28</td><td>28</td><td>28</td><td>28</td><td>28</td><td>28</td><td>28</td><td>28</td></tr><tr><td>Tip diameter mm</td><td>0.15</td><td>0.20</td><td>0.25</td><td>0.30</td><td>0.35</td><td>0.40</td><td>0.45</td><td>0.50</td><td>0.55</td><td>0.60</td><td>0.70</td><td>0.80</td></tr></table>		Tamaño	Diámetro	Color	Tamaño	Diámetro	Color	15	0.15	Blanco	45	0.45	Blanco	20	0.20	Amarillo	50	0.50	Amarillo	25	0.25	Rojo	55	0.55	Rojo	30	0.30	Azul	60	0.60	Azul	35	0.35	Verde	70	0.70	Verde	40	0.40	Negro	80	0.80	Negro	Standard specification	15	20	25	30	35	40	45	50	55	60	70	80	Length mm	28	28	28	28	28	28	28	28	28	28	28	28	Tip diameter mm	0.15	0.20	0.25	0.30	0.35	0.40	0.45	0.50	0.55	0.60	0.70	0.80
Tamaño	Diámetro	Color	Tamaño	Diámetro	Color																																																																														
15	0.15	Blanco	45	0.45	Blanco																																																																														
20	0.20	Amarillo	50	0.50	Amarillo																																																																														
25	0.25	Rojo	55	0.55	Rojo																																																																														
30	0.30	Azul	60	0.60	Azul																																																																														
35	0.35	Verde	70	0.70	Verde																																																																														
40	0.40	Negro	80	0.80	Negro																																																																														
Standard specification	15	20	25	30	35	40	45	50	55	60	70	80																																																																							
Length mm	28	28	28	28	28	28	28	28	28	28	28	28																																																																							
Tip diameter mm	0.15	0.20	0.25	0.30	0.35	0.40	0.45	0.50	0.55	0.60	0.70	0.80																																																																							
ENVASE	Contiene 200 unidades.																																																																																		
VIDA UTIL	N/A.																																																																																		
FINALIDAD PREVISTA	Sólo para uso dental. Se usa para secar el canal radicular.																																																																																		
MÉTODO DE ESTERILIZACIÓN	Producto estéril. No se garantiza la esterilidad después de haber abierto el envase. Método de esterilización: Rayos gamma. No esterilizar en autoclave o por calor.																																																																																		
ALMACENAMIENTO Y TRANSPORTE	No special storage and transport conditions required.																																																																																		

TECHNICAL DATA SHEET

BESTDENT PAPER TIPS																																																																																			
REFERENCES	See annex 1.																																																																																		
MANUFACTURER	Shanghai Xingyu Medical Equipment Co., Ltd. 3 Floor, Building 21, No. 3825, Xin Zhuan Road 201601 Shanghai, China																																																																																		
CLASSIFICATION	Medical device class IIa, according to Directive 93/42/EEC.																																																																																		
DESCRIPTION	Absorbent paper tips for drying root canals prepared with conventional instruments prior to obturation. ISO colour coding for easy identification.																																																																																		
SPECIFICATION	<table><tr><td>Tamaño</td><td>Diámetro</td><td>Color</td><td>Tamaño</td><td>Diámetro</td><td>Color</td></tr><tr><td>15</td><td>0.15</td><td>Blanco</td><td>45</td><td>0.45</td><td>Blanco</td></tr><tr><td>20</td><td>0.20</td><td>Amarillo</td><td>50</td><td>0.50</td><td>Amarillo</td></tr><tr><td>25</td><td>0.25</td><td>Rojo</td><td>55</td><td>0.55</td><td>Rojo</td></tr><tr><td>30</td><td>0.30</td><td>Azul</td><td>60</td><td>0.60</td><td>Azul</td></tr><tr><td>35</td><td>0.35</td><td>Verde</td><td>70</td><td>0.70</td><td>Verde</td></tr><tr><td>40</td><td>0.40</td><td>Negro</td><td>80</td><td>0.80</td><td>Negro</td></tr></table> Length: 28-+2mm. Diameter means the diameter of the tip of the product. (Taper 0.02, Taper 0.04, Taper 0.6) <table><tr><td>Standard specification</td><td>15</td><td>20</td><td>25</td><td>30</td><td>35</td><td>40</td><td>45</td><td>50</td><td>55</td><td>60</td><td>70</td><td>80</td></tr><tr><td>Length mm</td><td>28</td><td>28</td><td>28</td><td>28</td><td>28</td><td>28</td><td>28</td><td>28</td><td>28</td><td>28</td><td>28</td><td>28</td></tr><tr><td>Tip diameter mm</td><td>0.15</td><td>0.20</td><td>0.25</td><td>0.30</td><td>0.35</td><td>0.40</td><td>0.45</td><td>0.50</td><td>0.55</td><td>0.60</td><td>0.70</td><td>0.80</td></tr></table>		Tamaño	Diámetro	Color	Tamaño	Diámetro	Color	15	0.15	Blanco	45	0.45	Blanco	20	0.20	Amarillo	50	0.50	Amarillo	25	0.25	Rojo	55	0.55	Rojo	30	0.30	Azul	60	0.60	Azul	35	0.35	Verde	70	0.70	Verde	40	0.40	Negro	80	0.80	Negro	Standard specification	15	20	25	30	35	40	45	50	55	60	70	80	Length mm	28	28	28	28	28	28	28	28	28	28	28	28	Tip diameter mm	0.15	0.20	0.25	0.30	0.35	0.40	0.45	0.50	0.55	0.60	0.70	0.80
Tamaño	Diámetro	Color	Tamaño	Diámetro	Color																																																																														
15	0.15	Blanco	45	0.45	Blanco																																																																														
20	0.20	Amarillo	50	0.50	Amarillo																																																																														
25	0.25	Rojo	55	0.55	Rojo																																																																														
30	0.30	Azul	60	0.60	Azul																																																																														
35	0.35	Verde	70	0.70	Verde																																																																														
40	0.40	Negro	80	0.80	Negro																																																																														
Standard specification	15	20	25	30	35	40	45	50	55	60	70	80																																																																							
Length mm	28	28	28	28	28	28	28	28	28	28	28	28																																																																							
Tip diameter mm	0.15	0.20	0.25	0.30	0.35	0.40	0.45	0.50	0.55	0.60	0.70	0.80																																																																							
PACKAGING	Contains 200 units.																																																																																		
SHELF LIFE	N/A.																																																																																		
INTENDED USE	For dental use only. Used to dry the root canal.																																																																																		
STERILIZATION	Sterile product. Sterility is not guaranteed after opening. Sterilisation method: Gamma rays. Do not autoclave or heat sterilise.																																																																																		
STORAGE AND TRANSPORT	No special storage and transport conditions required.																																																																																		

ANEXO/ANNEX I:

Referencia	Descripción
100437	PUNTAS PAPEL BESTDENT N.15
100438	PUNTAS PAPEL BESTDENT N.20
100439	PUNTAS PAPEL BESTDENT N.25
100440	PUNTAS PAPEL BESTDENT N.30
100441	PUNTAS PAPEL BESTDENT N.35
100442	PUNTAS PAPEL BESTDENT N.40
100443	PUNTAS PAPEL BESTDENT N.45
100444	PUNTAS PAPEL SURTIDO 15-40 BESTDENT
100445	PUNTAS PAPEL BESTDENT N.50
100446	PUNTAS PAPEL BESTDENT N.55
100447	PUNTAS PAPEL BESTDENT N.60
100448	PUNTAS PAPEL BESTDENT N.70
100449	PUNTAS PAPEL BESTDENT N.80
100450	PUNTAS PAPEL SURTIDO 45-80 BESTDENT
100451	PUNTAS PAPEL CONICIDAD 04 N.20 BESTDENT
100452	PUNTAS PAPEL CONICIDAD 04 N.25 BESTDENT
100453	PUNTAS PAPEL CONICIDAD 04 N.30 BESTDENT
100454	PUNTAS PAPEL CONICIDAD 04 N.15-40 BESTDE
100455	PUNTAS PAPEL CONICIDAD 06 N.20 BESTDENT
100456	PUNTAS PAPEL CONICIDAD 06 N.25 BESTDENT
100457	PUNTAS PAPEL CONICIDAD 06 N.30 BESTDENT
100458	PUNTAS PAPEL CONICIDAD 06 N.15-40 BESTDE
100459	PUNTAS PAPEL TAPER F1 BESTDENT
100460	PUNTAS PAPEL TAPER F2 BESTDENT
100461	PUNTAS PAPEL TAPER F3 BESTDENT
100462	PUNTAS PAPEL TAPER SURT. F1-F2-F3 BD
100463	PUNTAS PAPEL TAPER SURT.F4-F5 BESTDENT
100468	PUNTAS PAPEL REC R25 (100u) BD
100469	PUNTAS PAPEL REC R40 (100u) BD
100470	PUNTAS PAPEL REC R50 (100u) BD
100471	PUNTAS PAPEL REC SURTIDA(100u) BD