

FICHA TÉCNICA

LIMAS K FLEXIBLES N.15-40 PROCLINIC		
REFERENCIA	Ver anexo 1	
FABRICANTE / RESPONSABLE	SHENZHEN SUPERLINE TECHNOLOGY CO., LTD. Room 302, Nanshan Zhiyuan A4, No.1001 Xueyuan Avenue, Changyuan community, Taoyuan Street, Nanshan District, Shenzhen, 518000 Guangdong, P.R. China	
CLASIFICACIÓN DEL PRODUCTO	Producto Sanitario clase IIa, según la Directiva 93/42/CEE relativa a productos sanitarios.	
DESCRIPCIÓN	- Las limas con una gran flexibilidad progresiva, que permite que el instrumento se adapte a la anatomía del canal radicular hasta el ápice. - Más flexible que las limas 'K' convencionales.	
COMPOSICIÓN	Parte activa: Acero inoxidable. Tope: Goma. Mango: Plástico.	
ENVASE	Envase que contiene 6 unidades.	
FINALIDAD PREVISTA	Indicadas para el tratamiento de canales curvos.	
LIMPIEZA Y ESTERILIZACIÓN	Los dispositivos se deben examinar, limpiar y esterilizar antes de cada uso. Usar detergentes y desinfectantes adecuados para instrumentos médicos o dentales. Desechar el producto si está dañado o disminuye su rendimiento.	
ALMACENAMIENTO Y TRANSPORTE	No requiere condiciones especiales de almacenamiento y transporte.	

TECHNICAL DATA SHEET

PROCLINIC FLEXIBLE K FILES N.15-40		
REFERENCES	See Annex 1.	
MANUFACTURER	SHENZHEN SUPERLINE TECHNOLOGY CO., LTD. Room 302, Nanshan Zhiyuan A4, No.1001 Xueyuan Avenue, Changyuan community, Taoyuan Street, Nanshan District, Shenzhen, 518000 Guangdong, P.R. China	
CLASSIFICATION	Medical device class IIa, according to Directive 93/42/EEC on medical devices.	
DESCRIPTION	- Files with high progressive flexibility, which allows the instrument to adapt to the anatomy of the root canal up to the apex. - More flexible than conventional 'K' files.	
COMPOSITION	Active part: Stainless steel. Stop: Rubber. Handle: Plastic.	
PACKAGING	Package containing 6 units.	
INTENDED USE	Indicated for the treatment of curved canals.	
CLEANING AND STERILIZATION	Devices should be examined, cleaned and sterilized before each use. Use detergents and disinfectants suitable for medical or dental instruments. Discard the product if it is damaged or its performance decreases.	
STORAGE AND TRANSPORT	It does not require special storage and transportation conditions.	

ANEXO / ANNEX I

REFERENCIA REFERENCES	DESCRIPCIÓN DESCRIPTION
17689	LIMAS K FLEX N.15 (25mm) PROC 6u.
17690	LIMAS K FLEX N.30 (25mm) PROC 6u.
17691	LIMAS K FLEX N.20 (25mm) PROC 6u.
17692	LIMAS K FLEX N.25 (25mm) PROC 6u.
17693	LIMAS K FLEX N.35 (25mm) PROC 6u.
17694	LIMAS K FLEX N 20 (21mm) PROC 6u.
17695	LIMAS K FLEX N.40 (25mm) PROC 6u.
17696	LIMAS K FLEX N. 15-40 (25mm) PROC 6u.
17697	LIMAS K FLEX N.15 (21mm) PROC 6u.
17698	LIMAS K FLEX N 25 (21mm) PROC 6u.
17699	LIMAS K FLEX N 30 (21mm) PROC 6u.
17700	LIMAS K FLEX N 35 (21mm) PROC 6u.
17701	LIMAS K FLEX N 40 (21mm) PROC 6u.
17702	LIMAS K FLEX N. 15-40 (21mm) PROC 6u.
17703	LIMAS K FLEX N.30 (31mm) PROC 6u.
17704	LIMAS K FLEX N.20 (31mm) PROC 6u.
17705	LIMAS K FLEX N.15 (31mm) PROC 6u.
17706	LIMAS K FLEX N.25 (31mm) PROC 6u.
17707	LIMAS K FLEX N.35 (31mm) PROC 6u.
17708	LIMAS K FLEX N.40 (31mm) PROC 6u.
17710	LIMAS K FLEX N. 15-40 (31mm) PROC 6u.