

FICHA TÉCNICA

FUNDAS PROTECTORAS PARA PLACAS DE FÓSFORO PROCLINIC		
REFERENCIAS	20497: Nº 1 49991: Nº 2 49992: Nº 0	
FABRICANTE / RESPONSABLE	J. Ripoll S.L Río Tiétar, 20 28110, Algete, Madrid España	
CLASIFICACIÓN	Producto No Sanitario (Clase consumo).	
DESCRIPCIÓN	Fundas protectoras para placas de fósforo que garantizan una óptima calidad de imagen. Aportan comodidad y seguridad al paciente evitando contaminación cruzada. Características: -Fáciles de manejar por su diseño ergonómico. -Diseño especial en la parte superior que simplifica la apertura de los sobres. -Compatibles con los distintos tamaños existentes de placas.	
ESPECIFICACIÓN	20497: Nº 1: (2 x 4 cm). 49991: Nº 2: (3 x 4 cm). 49992: Nº 0: (2 x 3 cm).	
COMPOSICIÓN	Polipropileno.	
PRESENTACIÓN	20497: Nº 1: 100 unidades por caja. 49991: Nº 2: 300 unidades por caja. 49992: Nº 0: 100 unidades por caja.	
FINALIDAD	Cobertura de las placas de fósforo utilizadas para las radiografías, evitando la contaminación cruzada entre pacientes.	
ALMACENAMIENTO Y TRANSPORTE	Mantener los botes cerrados y en un lugar seco y fresco.	

TECHNICAL DATA SHEET

BARRIER ENVELOPES PROCLINIC		
REFERENCES	20497: Nº 1 49991: Nº 2 49992: Nº 0	
MANUFACTURER	J. Ripoll S.L Río Tiétar, 20 28110, Algete, Madrid (SPAIN)	
CLASSIFICATION	Non Sanitary Product (Consumer Class).	
DESCRIPTION	Protective covers for phosphor plates that guarantee optimum image quality. They provide comfort and safety to the patient avoiding cross contamination. Features: -Easy to handle due to its ergonomic design. -Special design in the upper part that simplifies the opening of the envelopes. -Compatible with the different existing plate sizes.	
SPECIFICATION	20497: Nº 1: (2 x 4 cm). 49991: Nº 2: (3 x 4 cm). 49992: Nº 0: (2 x 3 cm).	
COMPOSITION	Polypropylene.	
PACKAGING	20497: Nº 1: 100 units per box. 49991: Nº 2: 300 units per box. 49992: Nº 0: 100 units per box.	
INTENDED USE	Coverage of the phosphor plates used for the X-rays, avoiding cross contamination between patients.	
HANDLING AND STORAGE CONDITIONS	Keep the cans closed and in a dry and cool place.	