

FICHA TÉCNICA

TIJERA IRIS		
REFERENCIAS	59833: Tijera Iris Curva Super Corte. 59834: Tijera Iris Recta Super Corte.	
FABRICANTE / RESPONSABLE	ASA DENTAL SPA Via Valenzana, 60 55054 MASSAROSA – BOZZANO LUCCA - ITALY (IT)	
CLASIFICACIÓN	Producto Sanitario Clase I, según la Directiva 93/42/CEE relativa a productos sanitarios.	
DESCRIPCIÓN	Tijera Iris para cirugía autoafilantes. Características: Ideal para realizar cortes y retirar suturas. Super corte.	
ESPECIFICACIÓN	Longitud: 11 cm tanto para las tijeras rectas como curvas.	
COMPOSICIÓN	Carburo de tungsteno y acero inoxidable.	
VIDA ÚTIL	La vida útil del instrumento depende de la frecuencia de uso, el cuidado del usuario y los métodos de reprocesamiento adecuados.	
PRESENTACIÓN	Contiene una unidad.	
FINALIDAD	Herramienta ideal para retirar suturas en zonas anteriores o realizar pequeños cortes de precisión.	
ESTERILIZACIÓN	Todos los instrumentos se deben examinar, limpiar y esterilizar antes de cada uso. Esterilizar a 130 ° en autoclave durante al menos 4 minutos, seguido de un secado adecuado.	
ALMACENAMIENTO Y TRANSPORTE	No requiere condiciones especiales de almacenamiento y transporte.	

TECHNICAL DATA SHEET

IRIS SIZE		
REFERENCE	59833: Curved Iris Scissors Super Cut. 59834: Super Cut Straight Iris Scissors.	
MANUFACTURER	ASA DENTAL SPA Via Valenzana, 60 55054 MASSAROSA – BOZZANO LUCCA - ITALY (IT)	
CLASSIFICATION	Class I Medical Device, Directive 93/42/EEC.	
DESCRIPTION	Iris scissors for self-sharpening surgery. Features: Ideal for cutting and suture removal. Super cut.	
SPECIFICATION	Length: 11 cm for both straight and curved scissors.	
COMPOSITION	Tungsten carbide and stainless steel.	
SHELF LIFE	Instrument lifetime depends on frequency of use, user care and proper reprocessing methods.	
PRESENTATION	1 unit.	
INTENDED USE	Ideal tool for removing sutures in anterior areas or making small precision cuts.	
STERILIZATION	All instruments should be examined, cleaned and sterilized before each use. Sterilize at 130° in autoclave for at least 4 minutes, followed by adequate drying.	
STORAGE AND TRANSPORT	It does not require special storage and transportation conditions.	