

FICHA TÉCNICA

SUTURA SEDA BESTDENT		
REFERENCIAS	69355: SUTURA SEDA 3/0 C-14. 69356: SUTURA SEDA 4/0 C-14. 69357: SUTURA SEDA 3/0 TC-12. 69358: SUTURA SEDA 3/0 TC-14. 69359: SUTURA SEDA 3/0 TC-16. 69360: SUTURA SEDA 4/0 TC-14. 69361: SUTURA SEDA 3/0 TB-15. 69362: SUTURA SEDA 4/0 TB-12. 69363: SUTURA SEDA 4/0 TB-15.	
FABRICANTE / RESPONSABLE	SMI AG Steinerberg 8, 4780 ST. VITH – BELGIUM Phone.: + 32 - (0)80 - 227292 Fax : + 32 – (0)80 - 226918 E-mail : info@sutures.be	
CLASIFICACIÓN	Producto Sanitario Clase IIb, de acuerdo con la Directiva 93/42/CEE, cumpliendo el reglamento 607/2023 (UE), relativo a los nuevos pedidos transitorios para la aplicación del reglamento 2017/745 (UE).	
DESCRIPCIÓN	- Sutura de seda natural trenzada Bestdent que proporciona un excelente y continuado soporte tisular hasta 3 meses. - Revestida de cera y silicona. - Alta fuerza tensil, no reabsorbible. - Suave con los tejidos blandos, produciendo el mínimo trauma posible. - Esterilizadas, en envase hermético para conservar las condiciones de esterilidad. - Color negro.	
ESPECIFICACIÓN	<u>Medidas:</u> 69355: (1/2-C 22MM) 69356: (1/2-C 22MM) 69357: (1/2-T 20MM) 69358: (1/2-T 22MM) 69359: (1/2-T 26MM) 69360: (1/2-T 22MM) 69361: (3/8-T 19MM) 69362: (3/8-T 16MM) 69363: (3/8-T 19MM)	
COMPOSICIÓN	Seda natural, cera y silicona.	
PRESENTACIÓN	12 unidades de 75 cm por caja.	
FINALIDAD PREVISTA	Suturar o ligar tejidos.	

FICHA TÉCNICA

INSTRUCCIONES DE USO	Ver instrucciones en el producto.
ALMACENAMIENTO Y TRANSPORTE	El producto debe protegerse de la acción directa de la luz solar y del calor y conservarse en el envase original en lugares limpios y secos, a una temperatura de entre 5 °C y 30 °C. No utilizar después de la fecha de caducidad.

TECHNICAL DATA SHEET

SILK SUTURE BESTDENT		
REFERENCES	69355: SILK SUTURE 3/0 C-14. 69356: SILK SUTURE 4/0 C-14. 69357: SILK SUTURE 3/0 TC-12. 69358: SILK SUTURE 3/0 TC-14. 69359: SILK SUTURE 3/0 TC-16. 69360: SILK SUTURE 4/0 TC-14. 69361: SILK SUTURE 3/0 TB-15. 69362: SILK SUTURE 4/0 TB-12. 69363: SILK SUTURE 4/0 TB-15	
MANUFACTURER	SMI AG Steinerberg 8, 4780 ST. VITH – BELGIUM Phone.: + 32 - (0)80 - 227292 Fax: + 32 – (0)80 - 226918 E-mail: info@sutures.be	
CLASSIFICATION	Class IIb Medical Device, in accordance with Directive 93/42/EEC, complying with regulation 607/2023 (EU), regarding new transitional orders for the application of regulation 2017/745 (EU).	
DESCRIPTION	- Bestdent braided natural silk suture that provides excellent and continuous tissue support for up to 3 months. - Wax and silicone coated. - High tensile strength, non-resorbable. - Gentle on soft tissues, producing the minimum possible trauma. - Sterilized, in hermetic packaging to preserve sterile conditions. - Black color.	
SPECIFICATION	<u>MEASURES:</u> 69355: (1/2-C 22MM) 69356: (1/2-C 22MM) 69357: (1/2-T 20MM) 69358: (1/2-T 22MM) 69359: (1/2-T 26MM) 69360: (1/2-T 22MM) 69361: (3/8-T 19MM) 69362: (3/8-T 16MM) 69363: (3/8-T 19MM)	
COMPOSITION	Natural silk, wax and silicone.	
PACKAGING	12 units of 75 cm per box	
INTENDED USE	Suture or ligate tissues.	

TECHNICAL DATA SHEET

INSTRUCTIONS FOR USE	See instructions on the product.
STORAGE AND TRANSPORT	The product should be protected from direct sunlight and heat and stored in the original container in a clean and dry place a temperature between 5 °C and 30 °C. Do not use after the expiration date.