

FICHA TÉCNICA

BLOCK DE MEZCLA PVC BESTDENT		
REFERENCIA	69675	
FABRICANTE / RESPONSABLE	LARIDENT SRL Via Lamaneigra, 12R, 16030 Ferrada GE, Italia	
CLASIFICACIÓN	Producto No sanitario (Clase consumo).	
DESCRIPCIÓN	- Blocks de mezcla de color transparente, con base antideslizante. - Medidas de 8X7 CM.	
COMPOSICIÓN	Realizados en PVC.	
ENVASE	1 unidad con 100 hojas.	
FINALIDAD PREVISTA	Bloques para realizar la mezcla de materiales dentales de impresión.	
INFORMACIÓN SOBRE LA SEGURIDAD	Producto desechable. Suministrado no estéril. No resiste a los hidrocarburos, resiste al alcohol y al aceite. Resistencia a la temperatura hasta 60 °C.	
CONDICIONES DE ALMACENAMIENTO Y TRANSPORTE	Conservar en lugar fresco y seco, alejado de la luz solar directa. Proteger de la humedad.	

TECHNICAL DATA SHEET

MIXING PAD PVC BESTDENT		
REFERENCES	69675	
MANUFACTURER / RESPONSIBLE	LARIDENT SRL Via Lamaneigra, 12R, 16030 Ferrada GE, Italia	
CLASSIFICATION	Non-sanitary product (consumer class).	
DESCRIPTION	- Blocks of transparent color mix, with non-slip base. - Measures 8X7 CM.	
COMPOSITION	Made of PVC.	
PACKAGE	1 unit with 100 sheets.	
INTENDED USE	Blocks for mixing dental impression materials.	
SAFETY INFORMATION	Disposable product. Supplied non-sterile. Not resistant to hydrocarbons, resistant to alcohol and oil. Temperature resistance up to 60 °C.	
STORAGE AND TRANSPORT CONDITIONS	Store in a cool, dry place away from direct sunlight. Protect from humidity.	