

FICHA TÉCNICA

POSTES ANATÓMICOS DE FIBRA DE VIDRIO PROCLINIC		
REFERENCIAS	7790	
FABRICANTE / RESPONSABLE	BIOLOREN S.R.L Via Alessandro Volta, 59 21047 Saronno Italy info@bioloren.com	
CLASIFICACIÓN	Producto Sanitario Clase IIa, Directiva 93/42/CEE.	
DESCRIPCIÓN	Postes anatómicos de fibra de vidrio translúcidos. Radiopacos. Características: Módulo de elasticidad similar a la dentina, ausencia de fractura radicular, estructura con fibra longitudinal, ausencia de oxidación y corrosión, forma cilindro-cónica. Ventajas: Fácil remoción del poste en caso de retratamiento, estabilidad y biocompatibilidad, ausencia de fracturas en el poste, uniforme distribución del estrés oclusal.	
COMPOSICIÓN	Fibra especial de vidrio y matriz resinosa epóxida. Composición en base a sulfato de bario.	
PRESENTACIÓN	KIT: 30 postes en 4 tamaños + 4 fresas.	
FINALIDAD	Ideales para las restauraciones estéticas de cerámica y el uso con cementos de polimerización dual. Radiopacos.	
INSTRUCCIONES DE USO	Ver instrucciones en el producto.	
ALMACENAMIENTO Y TRANSPORTE	No se requieren condiciones especiales de almacenamiento y transporte.	

TECHNICAL DATA SHEET

FIBERGLASS ANATOMICAL POSTS PROCLINIC		
REFERENCES	7790	
MANUFACTURER	BIOLOREN S.R.L Via Alessandro Volta, 59 21047 Saronno Italy info@bioloren.com	
CLASSIFICATION	Class IIa Medical Device, Directive 93/42/EEC.	
DESCRIPTION	Translucent fiberglass anatomical posts. Radiopaque. Features: Modulus of elasticity similar to dentin, absence of root fracture, structure with longitudinal fiber, absence of oxidation and corrosion, cylindrical-conical shape. Advantages: Easy removal of the post in case of retreatment, stability and biocompatibility, absence of fractures in the post, uniform distribution of occlusal stress.	
COMPOSITION	Special fiberglass and epoxy resin matrix. Composition based on barium sulfate.	
PACKAGING	KIT: 30 posts in 4 sizes + 4 drills.	
INTENDED USE	Ideal for esthetic ceramic restorations and use with dual-curing cements. Radiopaque.	
INSTRUCTIONS FOR USE	See instructions on the product.	
HANDLING AND STORAGE CONDITIONS	No special storage and transport conditions are required.	