

## FICHA TÉCNICA

|                                    |                                                                                                                                                                                                                                                                                                           |                                                                                    |
|------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------|
| <b>MECHERO DE ALCOHOL METÁLICO</b> |                                                                                                                                                                                                                                                                                                           |  |
| <b>REFERENCIA</b>                  | 80221                                                                                                                                                                                                                                                                                                     |                                                                                    |
| <b>FABRICANTE / RESPONSABLE</b>    | <b>CHL MEDICAL SOLUTIONS, Srl</b><br>Via Felice Cavallotti, 15<br>20122 Milán - Italia.                                                                                                                                                                                                                   |                                                                                    |
| <b>CLASIFICACIÓN</b>               | Producto No Sanitario (Clase consumo).                                                                                                                                                                                                                                                                    |                                                                                    |
| <b>DESCRIPCIÓN</b>                 | Mechero de alcohol metálico con mecha incluida.                                                                                                                                                                                                                                                           |                                                                                    |
| <b>ESPECIFICACIÓN</b>              | Diámetro: 7cm. Alto: 4 cm.                                                                                                                                                                                                                                                                                |                                                                                    |
| <b>COMPOSICIÓN</b>                 | Acero inoxidable.                                                                                                                                                                                                                                                                                         |                                                                                    |
| <b>ENVASE</b>                      | 1 unidad.                                                                                                                                                                                                                                                                                                 |                                                                                    |
| <b>FINALIDAD PREVISTA</b>          | <ul style="list-style-type: none"><li>- Se utiliza tanto en laboratorios como en clínicas dentales para realizar experimentos como calentar, esterilizar instrumentos, etc.</li><li>- Principalmente es usado para calentar la cera y poderla manipular en la realización de prótesis dentales.</li></ul> |                                                                                    |
| <b>ALMACENAMIENTO Y TRANSPORTE</b> | No requiere condiciones especiales de almacenamiento y transporte.                                                                                                                                                                                                                                        |                                                                                    |

## TECHNICAL DATA SHEET

|                         |                                                                                                                                                                                                                                                                                        |                                                                                     |
|-------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| ALCOHOL BURNER METALLIC |                                                                                                                                                                                                                                                                                        |  |
| REFERENCE               | 80221                                                                                                                                                                                                                                                                                  |                                                                                     |
| MANUFACTURER            | CHL MEDICAL SOLUTIONS, Srl<br>Via Felice Cavallotti, 15<br>20122 Milán - Italia.                                                                                                                                                                                                       |                                                                                     |
| CLASSIFICATION          | Non-Sanitary Product (Consumer Class).                                                                                                                                                                                                                                                 |                                                                                     |
| DESCRIPTION             | Metallic alcohol lighter with wick included.                                                                                                                                                                                                                                           |                                                                                     |
| SPECIFICATION           | Diameter: 7 cm. Height: 4 cm.                                                                                                                                                                                                                                                          |                                                                                     |
| COMPOSITION             | Stainless steel.                                                                                                                                                                                                                                                                       |                                                                                     |
| PACKAGE                 | 1 unit.                                                                                                                                                                                                                                                                                |                                                                                     |
| INTENDED USE            | <ul style="list-style-type: none"><li>- It is used both in laboratories and in dental clinics for experiments such as heating, sterilizing instruments, etc.</li><li>- It is mainly used to heat the wax to be able to manipulate it in the production of dental prostheses.</li></ul> |                                                                                     |
| STORAGE AND TRANSPORT   | It does not require special storage and transportation conditions.                                                                                                                                                                                                                     |                                                                                     |