

FICHA TÉCNICA

LÉNTULOS N.º 25-40 BESTDENT		
REFERENCIAS	87958: N.25 BD (25mm). 87959: N.30 BD (25mm). 87960: N.35 BD (25mm). 87961: N.40 BD (25mm).	
FABRICANTE / RESPONSABLE	Shenzhen Perfect Medical Instruments Co.,Ltd Room 103, Building 3, No. 2, Weiqun Road, 4th Community, Henggang Street, Longgang District, Shenzhen, 518115, Guangdong, P.R. China	
CLASIFICACIÓN	Producto Sanitario Clase IIa, de acuerdo con la Directiva 93/42/CEE, cumpliendo el reglamento 607/2023 (UE), relativo a los nuevos periodos transitorios para la aplicación del reglamento 2017/745 (UE).	
DESCRIPCIÓN	LÉNTULOS N.º 25-40. <u>Características:</u> - Su flexibilidad permite rellenar los conductos y fondo de la cavidad pulpar. - Esterilizables en autoclave.	
ESPECIFICACIONES	Medidas: 25 mm.	
ENVASE	Contiene 4 unidades.	
COMPOSICIÓN	Fabricadas en aleación de acero de gran calidad.	
FINALIDAD PREVISTA	Instrumento de endodoncia que facilita el sellado y cementación de conductos radiculares.	
ALMACENAMIENTO Y TRANSPORTE	No requiere condiciones especiales de almacenamiento y transporte.	

TECHNICAL DATA SHEET

LENTULOS N.º 25-40 BESTDENT		
REFERENCES	87958: N.25 BD (25mm). 87959: N.30 BD (25mm). 87960: N.35 BD (25mm). 87961: N.40 BD (25mm).	
MANUFACTURER / RESPONSIBLE	Shenzhen Perfect Medical Instruments Co.,Ltd Room 103, Building 3, No. 2, Weiqun Road, 4th Community, Henggang Street, Longgang District, Shenzhen, 518115, Guangdong, P.R. China	
CLASSIFICATION	Class IIa Medical Device, in accordance with Directive 93/42/EEC, complying with regulation 607/2023 (EU), regarding the new transitional periods for the application of regulation 2017/745 (EU).	
DESCRIPTION	LENTULOS N. º 25-40. <u>Features:</u> - Its flexibility allows filling the canals and bottom of the pulp cavity. - Sterilizable in autoclave.	
SPECIFICATIONS	Measurements: 25 mmm.	
PACKAGE	Contains 4 units.	
COMPOSITION	Made of high-quality steel alloy.	
INTENDED USE	Endodontic instrument that facilitates the sealing and cementation of root canals.	
STORAGE AND TRANSPORT	It does not require special storage and transportation conditions.	