

FICHA TÉCNICA

SET DISCOS DE PULIR PROCLINIC		
REFERENCIAS	98080: Ø 10 MM 98074: Ø 14 MM	
FABRICANTE / RESPONSABLE	EVE Ernst Vetter GmbH Neureutstr. 6- 75210 Keltern Germany Tel.: +49 7231 97 77 -0 Fax: +49 7231 97 77 -99 E-Mail: info@eve-rotary.com	
CLASIFICACIÓN	Producto Sanitario Clase IIa, conforme el Reglamento (UE) 2017/745.	
DESCRIPCIÓN	Discos de plástico para un excelente acabado y pulido. Muy flexibles, recubiertos con óxido de aluminio. Disponibles en cuatro granos de abrasividad y dos diámetros (10 mm y 14 mm) para cubrir todas las necesidades clínicas.	
ESPECIFICACIÓN	<u>Tamaño del grano:</u> Grano grueso: 50-100 micras. Grano medio: 30-50 micras. Grano fino: 15-30 micras. Grano extrafino: 2-15 micras.	
COMPOSICIÓN	Recubierto por un lado con óxido de aluminio.	
PRESENTACIÓN	CAJA 100 unidades (25 unidades x 4 granos) + 1 mandril.	
FINALIDAD	Especialmente diseñados para el pulido interdental, el sistema de Discos de Pulido Proclinic ofrece una solución completa para obtener un brillo superior en cada restauración.	
INSTRUCCIONES DE USO	Ver instrucciones en el producto.	
ALMACENAMIENTO Y TRANSPORTE	No requiere condiciones especiales de almacenamiento y transporte.	

TECHNICAL DATA SHEET

POLISHING DISC SET PROCLINIC		
REFERENCES	98080: Ø 10 MM 98074: Ø 14 MM	
MANUFACTURER	EVE Ernst Vetter GmbH Neureutstr. 6- 75210 Keltern Germany Tel.: +49 7231 97 77 -0 Fax: +49 7231 97 77 -99 E-Mail: info@eve-rotary.com	
CLASSIFICATION	Class IIa Medical Device, in accordance with Regulation (EU) 2017/745.	
DESCRIPTION	Plastic discs for an excellent finish and polish. Very flexible, coated with aluminum oxide. Available in four abrasiveness grains and two diameters (10mm and 14mm) to cover all clinical needs.	
SPECIFICATION	<u>Grain size:</u> Coarse grain: 50-100 microns. Medium grain: 30-50 microns. Fine grain: 15-30 microns. Extra fine grain: 2-15 microns.	
COMPOSITION	Coated on one side with aluminium oxide	
PACKAGING	BOX 100 units (25 units x 4 grains) + 1 mandrel.	
INTENDED USE	Specially designed for interdental polishing, the Proclinic Polishing Disc system offers a complete solution to obtain a superior shine in every restoration.	
INSTRUCTIONS FOR USE	See instructions on the product.	
HANDLING AND STORAGE CONDITIONS	It does not require special storage and transport conditions.	