

FICHA TÉCNICA

PASTA DE PULIR PARA ACRILICO PROCLINIC		
REFERENCIAS	H02162	
FABRICANTE / RESPONSABLE	MEGADENTA Dentalprodukte GmbH Carl-Eschebach-Strasse 1A D-01454 Radeberg Alemania	
CLASIFICACIÓN	Producto No Sanitario (Clase Consumo).	
DESCRIPCIÓN	Pasta de Pulir para acrílico.	
ENVASE	Bote de 500 gr.	
FINALIDAD PREVISTA	Pasta de pulido Proclinic para abrillantar/pulir resinas acrílicas y metal.	
INSTRUCCIONES DE USO	Pre-pulido como de costumbre, prótesis con el producto con brillo y lana pulida o similar. Pulido. Luego enjuague bajo el agua. Agitar antes de usar	
ALMACENAMIENTO Y TRANSPORTE	No requiere condiciones especiales de almacenamiento y transporte.	

TECHNICAL DATA SHEET

POLISHING LIQUID FOR ACRYLIC AND METALS PROCLINIC			
REFERENCES	H02162		
MANUFACTURER / RESPONSIBLE	MEGADENTA Dentalprodukte GmbH Carl-Eschebach-Strasse 1A D-01454 Radeberg Germany		
CLASSIFICATION	Non-Sanitary Product (Consumer Class).		
DESCRIPTION	Polishing paste for acrylic.		
PACKAGE	Jar of 500 gr.		
INTENDED USE	Proclinic polishing paste for polishing/polishing acrylic resins and metal.		
INSTRUCCIONS FOR USE	Pre polishing as usual, denture with moistening this product and using woolen bobbing or similar polishing. Subsequently rinse under water		
STORAGE AND TRANSPORT	It does not require special storage and transportation conditions.		