

FICHA TÉCNICA

FRESA TUNGSTENO PROCLINIC		FOTOGRAFIA VER ANEXO 2
REFERENCIAS	Ver anexo 1.	
FABRICANTE / RESPONSABLE	Dedeco International, Inc. 11617 State Route 97 Long Eddy, NY 12760 Telephone No.: 845-887-4840	
CLASIFICACIÓN	Producto no sanitario.	
DESCRIPCIÓN	Fresa de tungsteno.	
COMPOSICIÓN	Carburo de tungsteno, cobalto, carburo de tantato, carburo de cromo, carburo de titanio, carburo de niobio.	
PRESENTACIÓN	1 unidad por envase.	
ESTERILIZACIÓN	Todos los instrumentos se deben examinar, limpiar y esterilizar antes de cada uso. Usar detergentes y desinfectantes adecuados para instrumentos médicos o dentales. El final de la vida útil de los instrumentos viene determinado por el desgaste y los daños durante el uso, y los instrumentos deben inspeccionarse para detectar defectos durante el proceso de limpieza preliminar. Esterilizable en autoclave o calor seco.	
FINALIDAD PREVISTA	Herramientas de corte utilizadas para cortar, taladrar, esmerilar y terminar superficies.	
ALMACENAMIENTO Y TRANSPORTE	No requiere condiciones especiales de almacenamiento y transporte.	

TECHNICAL DATA SHEET

CARBIDE BURS PROCLINIC		PHOTOGRAPHY SEE ANEX 2
REFERENCES	See anex 1.	
MANUFACTURER	Dedeco International, Inc. 11617 State Route 97 Long Eddy, NY 12760 Telephone No.: 845-887-4840	
CLASSIFICATION	Non-medical device.	
DESCRIPTION	Tungsten burs.	
COMPOSITION	Tungsten Carbide, Cobalt, Tantatum Carbide, Chromium Carbide, Titanium Carbide, Niobium Carbide.	
PACKAGING	1 unit per package.	
STERILIZATION	All instruments should be examined, cleaned and sterilized before each use. Use detergents and disinfectants suitable for medical or dental instruments. The end of life of instruments is determined by wear and damage during use, and instruments should be inspected for defects during the preliminary cleaning process. Autoclavable or dry heat sterilizable.	
INTENDED USE	Cutting tools used for cutting, drilling, grinding and finishing surfaces.	
HANDLING AND STORAGE CONDITIONS	It does not require special storage and transport conditions.	




ANEXO 1 / ANNEX 1:

REFERENCIAS / REFERENCES:

H20710	FRESA TUNGSTENO PM 71E-023
H20714	FRESA TUNGSTENO PM 79E-040
H20720	FRESA TUNGSTENO PM 137E-023
H20734	FRESA TUNGSTENO PM 351E-060
H20736	FRESA TUNGSTENO PM 79EF-040
H20742	FRESA TUNGSTENO PM 251EF-060
H20746	FRESA TUNGSTENO PM 257EF-023
H20754	FRESA TUNGSTENO PM 79FST-040
H20756	FRESA TUNGSTENO PM 129FST-023
H20758	FRESA TUNGSTENO PM 136FST-016
H20760	FRESA TUNGSTENO PM 138FST-023
H20762	FRESA TUNGSTENO PM 139FST-023
H20764	FRESA TUNGSTENO PM 251XG-060
H20770	FRESA TUNGSTENO PM 251GSQ-060
H20772	FRESA TUNGSTENO PM 79SG-045
H20774	FRESA TUNGSTENO PM 79SG-060
H20712	FRESA TUNGSTENO PM 77E-029
H20716	FRESA TUNGSTENO PM 79E-045
H20718	FRESA TUNGSTENO PM 79E-060
H20722	FRESA TUNGSTENO PM 138E-016
H20724	FRESA TUNGSTENO PM 251E-060L
H20726	FRESA TUNGSTENO PM 257E-023
H20728	FRESA TUNGSTENO PM 261E-023
H20730	FRESA TUNGSTENO PM 295E-023
H20732	FRESA TUNGSTENO PM 351E-040
H20768	FRESA TUNGSTENO PM 79SF-045
H20776	FRESA TUNGSTENO PM 351SG-060
H20766	FRESA TUNGSTENO PM 351XG-070
H20738	FRESA TUNGSTENO PM 79EF-045
H20740	FRESA TUNGSTENO PM 138EF-016
H20744	FRESA TUNGSTENO PM 251EF-060L
H20748	FRESA TUNGSTENO PM 261EF-023
H20750	FRESA TUNGSTENO PM 79SPF-040
H20751	FRESA TUNGSTENO PM 129SPF-023
H20752	FRESA TUNGSTENO PM 251SPF-060
H20778	FRESA TUNGSTENO PM 219-023

ANEXO 2 / ANNEX 2:

FOTOGRAFIAS / PHOTOGRAPHS:

H20710		H20758		H20718		H20766	
H20714		H20760		H20722		H20738	
H20720		H20762		H20724		H20740	
H20734		H20764		H20726		H20744	
H20736		H20770		H20728		H20748	
H20742		H20772		H20730		H20750	
H20746		H20774		H20732		H20751	
H20754		H20712		H20768		H20752	
H20756		H20716		H20776		H20778	