

FICHA TÉCNICA

PULIDOR CERAMICA		
REFERENCIAS	H211128: PULIDOR CERAMICA BLANCO 22X3MM CX10. H211129: PULIDOR CERAMICA BLANCO 22X3MM CX100. H211130: PULIDOR CERAMICA ROSA 22X3MM CX10. H211131: PULIDOR CERAMICA ROSA 22X3MM CX100.	
FABRICANTE / RESPONSABLE	EVE Ernst Vetter GmbH Neureutstr. 6- 75210 Keltern Germany Tel.: +49 7231 97 77 -0 Fax: +49 7231 97 77 -99 E-Mail: info@eve-rotary.com	
CLASIFICACIÓN	Producto No Sanitario (Clase protésico).	
DESCRIPCIÓN	- Para el prepulido para las cerámicas dentales se adapta perfectamente a la superficie a tratar. - Repasado y prepulido. - Para un mayor rendimiento del pulidor se recomienda trabajar entre 7.000-15.000 RPM/min.	
ESPECIFICACIÓN	H211128: 22X3MM CX10. H211129: 22X3MM CX100. H211130: 22X3MM CX10. H211131: 22X3MM CX100.	
PRESENTACIÓN	H211128 y H211130: 10 uds. H211129 y H211131: 100 uds.	
FINALIDAD	Pulidores en dos pasos para el repasado y pulido.	
ALMACENAMIENTO Y TRANSPORTE	No requiere condiciones especiales de almacenamiento y transporte.	

TECHNICAL DATA SHEET

CERAMIC POLISHER		
REFERENCES	H211128: CERAMIC POLISHER WHITE 22X3MM CX10. H211129: CERAMIC POLISH WHITE 22X3MM CX100. H211130: CERAMIC POLISH PINK 22X3MM CX10. H211131: CERAMIC POLISH PINK 22X3MM CX100.	
MANUFACTURER	EVE Ernst Vetter GmbH Neureutstr. 6- 75210 Keltern Germany Tel.: +49 7231 97 77 -0 Fax: +49 7231 97 77 -99 E-Mail: info@eve-rotary.com	
CLASSIFICATION	Non-Health Product (Prosthetic Class).	
DESCRIPTION	- For the pre-polishing of dental ceramics, it is perfectly adapted to the surface to be treated. - Reworking and pre-polishing. - For a better performance of the polisher it is recommended to work between 7.000-15.000 RPM/min.	
SPECIFICATION	H211128: 22X3MM CX10. H211129: 22X3MM CX100. H211130: 22X3MM CX10. H211131: 22X3MM CX100.	
PACKAGING	H211128 y H211130: 10 unit. H211129 y H211131: 100 unit.	
INTENDED USE	Two-step polishers for grinding and polishing.	
HANDLING AND STORAGE CONDITIONS	It does not require special storage and transport conditions.	