



FICHA TÉCNICA

PULIDOR CERAMICA DIAMANTE 26X2MM		
REFERENCIAS	H211132: PULIDOR CERAMICA DIAMANTE 26X2 GRUESO 1U. H211133: PULIDOR CERAMICA DIAMANTE 26X2 MEDIO 1U. H211134: PULIDOR CERAMICA DIAMANTE 26X2 FINO 1U.	
FABRICANTE / RESPONSABLE	EVE Ernst Vetter GmbH Neureutstr. 6- 75210 Keltern Germany Tel.: +49 7231 97 77 -0 Fax: +49 7231 97 77 -99 E-Mail: info@eve-rotary.com	
CLASIFICACIÓN	Producto No Sanitario (Clase protésico).	
DESCRIPCIÓN	- Sistema de pulido en 3 pasos (Grueso/Medio/Fino) para todos los materiales cerámicos. - Con diamante para conseguir la máxima eficiencia. - Antes del pulido realizar un repasado con el abrasivo con diamante. - Para un mayor rendimiento del pulidor se recomienda trabajar entre 7.000-12.000 RPM/min.	
ESPECIFICACIÓN	H211132: 26X2 GRUESO 1U. H211133: 26X2 MEDIO 1U. H211134: 26X2 FINO 1U.	
PRESENTACIÓN	1 unidad.	
FINALIDAD	Pulidores de diamante de alto brillo para todo tipo de materiales cerámicos.	
ALMACENAMIENTO Y TRANSPORTE	No requiere condiciones especiales de almacenamiento y transporte.	

TECHNICAL DATA SHEET

DIAMOND CERAMIC POLISHER 26X2MM		 H211132H211133H211134
REFERENCES	H211132: DIAMOND CERAMIC POLISHER 26X2 COARSE 1U. H211133: DIAMOND CERAMIC POLISHER 26X2 MEDIUM 1U. H211134: DIAMOND CERAMIC POLISHER 26X2 FINE 1U.	
MANUFACTURER	EVE Ernst Vetter GmbH Neureutstr. 6- 75210 Keltern Germany Tel.: +49 7231 97 77 -0 Fax: +49 7231 97 77 -99 E-Mail: info@eve-rotary.com	
CLASSIFICATION	Non-Health Product (Prosthetic Class).	
DESCRIPTION	- 3-step polishing system (Coarse/Medium/Fine) for all ceramic materials. - With diamond for maximum efficiency. - Prior to polishing, a diamond abrasive should be used. - For a better performance of the polisher it is recommended to work between 7.000-12.000 RPM/min.	
SPECIFICATION	H211132: 26X2 COARSE 1U. H211133: 26X2 MEDIUM 1U. H211134: 26X2 FINE 1U.	
PACKAGING	1 unit.	
INTENDED USE	High gloss diamond polishers for all types of ceramic materials.	
HANDLING AND STORAGE CONDITIONS	It does not require special storage and transport conditions.	