



FICHA TÉCNICA

PULIDOR DE CIRCONIO		
REFERENCIAS	H211141: PULIDOR CIRCONIO 26X2 MEDIO 1U. H211142: PULIDOR CIRCONIO 4X13 MEDIO1U. H211143: PULIDOR CIRCONIO 11X2 MEDIO 1U H211144: PULIDOR CIRCONIO 26X2 FINO 1U. H211145: PULIDOR CIRCONIO 4X13 FINO 1U. H211146: PULIDOR CIRCONIO 11X2 FINO 1U.	
FABRICANTE / RESPONSABLE	EVE Ernst Vetter GmbH Neureutstr. 6- 75210 Keltern Germany Tel.: +49 7231 97 77 -0 Fax: +49 7231 97 77 -99 E-Mail: info@eve-rotary.com	
CLASIFICACIÓN	Producto No Sanitario (Clase protésico).	
DESCRIPCIÓN	- Pulidor indicado para el pulido de restauraciones en circonio a volumen total. - Pulidor con alta concentración de diamante para conseguir una mayor eficiencia. - Antes del pulido se recomienda repasar con el abrasivo para circonio. - Para un mayor rendimiento del pulidor se recomienda trabajar entre 7.000-12.000 RPM/min.	
ESPECIFICACIÓN	H211141: 26X2 MEDIO 1U. H211142: 4X13 MEDIO1U. H211143: 11X2 MEDIO 1U. H211144: 26X2 FINO 1U. H211145: 4X13 FINO 1U. H211146: 11X2 FINO 1U.	
PRESENTACIÓN	1 unidad.	
FINALIDAD	Especialmente indicado para circonio monolítico.	
ALMACENAMIENTO Y TRANSPORTE	No requiere condiciones especiales de almacenamiento y transporte.	

TECHNICAL DATA SHEET

ZIRCONIA POLISHER			
REFERENCES	H211141: PULIDOR CIRCONIO 26X2 MEDIO 1U. H211142: PULIDOR CIRCONIO 4X13 MEDIO1U. H211143: PULIDOR CIRCONIO 11X2 MEDIO 1U H211144: PULIDOR CIRCONIO 26X2 FINO 1U. H211145: PULIDOR CIRCONIO 4X13 FINO 1U. H211146: PULIDOR CIRCONIO 11X2 FINO 1U.		
MANUFACTURER	EVE Ernst Vetter GmbH Neureutstr. 6- 75210 Keltern Germany Tel.: +49 7231 97 77 -0 Fax: +49 7231 97 77 -99 E-Mail: info@eve-rotary.com		
CLASSIFICATION	Non-Health Product (Prosthetic Class).		
DESCRIPTION	- Polisher indicated for the polishing of full volume zirconia restorations. - Polisher with high concentration of diamond for greater efficiency. - Before polishing it is recommended to polish with zirconia abrasive. - For a better performance of the polisher it is recommended to work between 7.000-12.000 RPM/min.		
SPECIFICATION	H211141: 26X2 MEDIO 1U. H211142: 4X13 MEDIO1U. H211143: 11X2 MEDIO 1U. H211144: 26X2 FINO 1U. H211145: 4X13 FINO 1U. H211146: 11X2 FINO 1U.		
PACKAGING	1 unit.		
INTENDED USE	Especially suitable for monolithic zirconia.		
HANDLING AND STORAGE CONDITIONS	It does not require special storage and transport conditions.		