

FICHA TÉCNICA

PULIDOR PARA ACRILICO		
REFERENCIAS	H211147: PULIDOR PARA ACRILICO 10X20 GRUESO CX10. H211148: PULIDOR PARA ACRILICO 10X20 MEDIO CX10. H211149: PULIDOR PARA ACRILICO 10X20 FINO CX10. H211150: PULIDOR PARA ACRILICO 10X24 GRUESO CX10. H211151: PULIDOR PARA ACRILICO 10X24 MEDIO CX10. H211152: PULIDOR PARA ACRILICO 10X24 FINO CX10.	
FABRICANTE / RESPONSABLE	EVE Ernst Vetter GmbH Neureutstr. 6- 75210 Keltern Germany Tel.: +49 7231 97 77 -0 Fax: +49 7231 97 77 -99 E-Mail: info@eve-rotary.com	
CLASIFICACIÓN	Producto No Sanitario (Clase protésico).	
DESCRIPCIÓN	- Pulidor con carburo de silicio para prepulir superficies de materiales acrílicos. - Para un mayor rendimiento del pulidor se recomienda trabajar entre 5.000-7.000 RPM/min.	
ESPECIFICACIÓN	H211147: 10X20 GRUESO CX10. H211148: 10X20 MEDIO CX10. H211149: 10X20 FINO CX10. H211150: 10X24 GRUESO CX10. H211151: 10X24 MEDIO CX10. H211152: 10X24 FINO CX10.	
COMPOSICIÓN	Pulidor con carburo de silicio.	
PRESENTACIÓN	10 unidades.	
FINALIDAD	Para alisado con el grano grueso. Pre-pulido con el grano medio y Pulido brillante con el grano fino.	
ALMACENAMIENTO Y TRANSPORTE	No requiere condiciones especiales de almacenamiento y transporte.	

TECHNICAL DATA SHEET

ACRYLIC POLISH			
REFERENCES	H211147: POLISHER FOR ACRYLIC 10X20 THICK CX10. H211148: POLISHER FOR ACRYLIC 10X20 MEDIUM CX10. H211149: POLISHER FOR ACRYLIC 10X20 FINE CX10. H211150: POLISHER FOR ACRYLIC 10X24 THICK CX10. H211151: POLISHER FOR ACRYLIC 10X24 MEDIUM CX10. H211152: POLISHER FOR ACRYLIC 10X24 FINE CX10.		
MANUFACTURER	EVE Ernst Vetter GmbH Neureutstr. 6- 75210 Keltern Germany Tel.: +49 7231 97 77 -0 Fax: +49 7231 97 77 -99 E-Mail: info@eve-rotary.com		
CLASSIFICATION	Non-Health Product (Prosthetic Class).		
DESCRIPTION	- Silicon carbide polisher for pre-polishing surfaces of acrylic materials. - For best performance of the polisher it is recommended to work between 5,000-7,000 RPM/min.		
SPECIFICATION	H211147: 10X20 THICK CX10. H211148: 10X20 MEDIUM CX10. H211149: 10X20 FINE CX10. H211150: 10X24 THICK CX10. H211151: 10X24 MEDIUM CX10. H211152: 10X24 FINE CX10.		
COMPOSITION	Silicon carbide polisher.		
PACKAGING	10 units.		
INTENDED USE	For smoothing with coarse grain. Pre-polishing with medium grain and bright polishing with fine grain.		
HANDLING AND STORAGE CONDITIONS	It does not require special storage and transport conditions.		