



**BEYOND USA Headquarters**  
12503 Exchange Dr. Suite 530  
Stafford, TX 77477 USA  
TEL +1-281-277-4352  
EMAIL: info@beyonddent.com

**BEYOND European Headquarters**  
Germany  
TEL +48-600-288-309  
EMAIL: europe@beyonddent.com

**BEYOND Asia Pacific Headquarters**  
Chaoyang District, Beijing  
China 100029  
EMAIL: asia@beyonddent.com



[www.beyonddent.com](http://www.beyonddent.com)





— Global Leader of In-Office Whitening

BEYOND DENTAL PRODUCTS  
SHOWCASE







# GLOBAL LEADER of In-Office Teeth Whitening!

Since our inception in Sunnyvale, California BEYOND has focused on bringing comfortable, long lasting whitening results to smiles across the globe. Over 40,000 systems worldwide attest to our impact since 2003.



**248,000+**

Satisfied Dentists



**27,000+**

Worldwide Treatment  
Locations



**2,500,000+**

Optimized Smiles





BEYOND is proud to announce that the Polus Advanced accelerator has garnered the top whitening system for 2019. This will be an unbroken chain of awards since the Polus series was first introduced in 2011.

**POLUS**<sup>®</sup>  
Advanced Ultra  
Whitening Accelerator



Top Whitening System Award

9th Year in a row!



Unique Technology of  
**Ultrasound + LightBridge™**

32<sup>ND</sup> ANNIVERSARY

**REALITY**

The information source for esthetic dentistry

2018

FirstLook

"The Polus Advanced Ultra light is definitely the star of this show, combining halogen, LED, and ultrasonic technologies."





## Polus<sup>®</sup> Advanced Ultra Whitening Accelerator


The BEYOND<sup>®</sup> Polus<sup>®</sup> Advanced Ultra is a feature-expanded version of our award-winning Polus series of whitening accelerators. In addition to the features that BEYOND customers have grown to know and love, the Advanced Ultra has an augmented ultrasound function that lets the practitioner reduce the treatment time while maintaining the same great results.

Our powerful halogen light source is filtered by two sets of 30 optical lenses, then transported via over 200,000 optical fibers to add with our LED light source. As this light accelerates the chemical breakdown of the whitening gel, the patent-pending ultrasound function utilizes a cavitation process that agitates and permits the whitening gel to target the enamel more comprehensively and in a shorter time frame. This means less chair time for the patient, less potential for sensitivity and the same great results you've come to expect with BEYOND.

### *Innovative and Power-Packed System*

- Upgraded light filtration system containing over 200,000 optical fibers, to filter out all harmful heat and ultraviolet (UV) light
- Combines high-intensity halogen lamps, medical LED lamps and ultrasonic generators to form a composite activation source of hydrogen peroxide
- Dual-arch power whitening option can be set to High, Medium, or Low light output settings for a tailored treatment that best suits each individual patient
- Single-tooth whitening option included for precision, tailored treatment results
- Smart Remote Control that allows patient awareness of procedure time
- Built in purifier for a cleaner and healthier environment
- Built-in instruction mode, that allows for a customized whitening procedure based on the client's needs
- Built-in, powerful LED curing light features a dual wavelength output to cure any dental composite on the market
- State of the art, full colored, touch screen with newly added multilingual feature
- Advanced temperature monitor and cooling system to prevent system from overheating





LightBridge™ Technology



Ultrasonic Technology



Touch User Interface (TUI)  
Control Panel



**VOICE INSTRUCTIONS**

Remote Controlled Voice  
Instruction Function



LED Curing Light /  
Spot Whitening Device



**POLUS**®  
Advanced Ultra  
Whitening Accelerator

Unprecedented

**Innovation and Flexibility**  
in Whitening Technology



before after

**A Brilliant Smile**  
in just **30 min.**

**Top Whitening System**  
in the world!

**DA**  
DENTAL ADVISOR  
2019  
TOP AWARD  
WINNER  
DENTAL ADVISOR  
TOP AWARD  
IN-OFFICE  
WHITENING

**DA**  
DENTAL ADVISOR  
2018  
TOP AWARD  
WINNER  
DENTAL ADVISOR  
TOP AWARD  
IN-OFFICE  
WHITENING

DENTAL ADVISOR  
TOP AWARD  
2014  
DENTAL ADVISOR  
TOP AWARD  
2013  
DENTAL ADVISOR  
TOP AWARD  
2012  
DENTAL ADVISOR  
TOP AWARD  
2011  
DENTAL ADVISOR  
TOP AWARD  
2017  
DENTAL ADVISOR  
TOP AWARD  
2016  
DENTAL ADVISOR  
TOP AWARD  
2015

**beyond**  
**POLUS**  
Advanced  
www.beyonddent.com

**A Brilliant Smile**  
in just **30 min.**

**beyond** International Inc.  
DENTAL ADVISOR TOP AWARD WINNER  
Products of USA

Wall Poster

Standing Poster

**beyond**  
**POLUS**  
Advanced  
Whitening Accelerator

**Innovative and Power-Packed System**  
Multifunctional system for greater treatment flexibility with LightBridge™ powered whitening accelerator (dual-arch and single-tooth whitening options), new ultrasonic technology, and built-in LED curing light.

**USER MANUAL**  
Instructions for Use

Instruction Manual

**A Brilliant Smile**  
in just **30 min.**

**beyond** International Inc.  
DENTAL ADVISOR TOP AWARD WINNER  
Products of USA

4-Fold Flyer

**A Brilliant Smile**  
in just **30 min.**

**beyond** International Inc.  
DENTAL ADVISOR TOP AWARD WINNER  
Products of USA

Wall Poster

Every BEYOND® Whitening Accelerator Comes With  
A Complete Marketing Package



BEYOND® Power Whitening is recognized as the  
**Safest and Most Effective** way to whiten your teeth



Before

After

Healthy,  
shiny results



Lasting color

Comfort in treatment



## Ultra

The **ultrasound** function utilizes a cavitation process that permits the whitening gel to target the enamel more comprehensively and in a shorter time frame. This means less chair time for the patient, less potential for sensitivity and the same great results you've come to expect with BEYOND®.

## Sleek

The BEYOND II Ultra has a **small footprint** and can be easily moved from operator to operator.



Whitening Made

Simple and

**beyond**® II Ultra Whitening Accelerator





## Simultaneous

Allows for the effective treatment of two patients at the same time.

**Second LED sold separately**

## Ionizer

A passive **air ionizer** gently removes odors from the room while unit is on.



# Affordable!





## **beyond**<sup>®</sup> II Ultra Whitening Accelerator

The chic and modern design of the BEYOND<sup>®</sup> II Whitening Accelerator works for any dental office, spa, or salon. The BEYOND II Ultra is the latest improvement to the popular BEYOND II whitening accelerator series. Along with the convenient features of simultaneous treatment capability, intuitive touch user interface and passive air ionizer, the LED light features our patent-pending ultrasound feature to aid in brevity of treatment time and little to no sensitivity.

### Professional in-office Whitening Treatment



Digital



Purifies



Adjustable



Powerful



**Elegant** – Chic and modern design that will be easily embraced by any dental practice.

**Adjustable** – Whitening option can be set to gentle, normal, or intense light output settings for a tailored treatment that best suits each individual patient.

**Digital** – Simple and quick patient preparation and system operation with digital display control board.


**Purifies** – Using negative ions, the built in air purifier reduces airborne pollutants that may exist within the practice, actively contributing to a cleaner and healthier environment.

When used in combination with BEYOND® proprietary whitening formula which contains hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) and silicon dioxide (SiO<sub>2</sub>) that act as a catalyst, along with other special ingredients, the end result is not only whiter teeth but shinier teeth with reduced prep and chair time.

- The BEYOND® II is both lightweight and portable for both travel and ease of use for room to room
- Whitening options can be set to Low, Medium and High light output settings for a tailored treatment that best suits each individual client
- Quick and simple client preparation and system operation with digital display control board
- Utilizing negative ions, the built-in air purifier reduces airborne pollutants that may exist within the treatment space, actively contributing to a cleaner and healthier environment




Wall Poster



Wall Poster



4-Fold Flyer




Instruction Manual

Every BEYOND® II Whitening Accelerator Comes With  
A Complete Marketing Package



***beyond***<sup>®</sup>  
**COMMAND**

# The First Generation of BEYOND<sup>®</sup> Whitening Accelerator



THE DENTAL  
ADVISOR  
++++ ½



12

## BEYOND<sup>®</sup> COMMAND<sup>®</sup> Whitening Accelerator

Stronger light-output than our original accelerator: Powerful halogen light for superior whitening results.

Advanced light-filtering system with over 13,500 optical fibers and 30 layers of coated optical glass removes harmful ultraviolet (UV) light and heat.

Tailored for use with our most popular whitening kits: BEYOND<sup>®</sup> II Complete Professional Whitening Kit, BEYOND<sup>®</sup> II Five Patient Economy Kit, and BEYOND<sup>®</sup> Blast. \*

\* Look for the COMMAND<sup>®</sup> sticker on compatible whitening kits





Brilliant smile, colorful life!



# Rebuild, Strengthen, and Whiten Formula

The BEYOND® Max5® Solution Treatment Kit is a revolutionary new product, combining our most popular in-office teeth whitening kit and a specially formulated RS gel to reduce known sensitivity of some patients. This kit includes our same dual-barrel syringe W Gel 1 with 35% hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) and our specially formulated RS Gel 2 that contains desensitizing ingredients. This combination of potassium nitrate and hydroxyapatite offers enamel regeneration and restoration. When combined, the strength of both gels reduce the hydrogen peroxide concentration to 6% and helps aid in the reduction of sensitivity, as well as remineralization of the enamel. This kit includes the essentials for chairside teeth bleaching of five patients. A typical treatment is three consecutive cycles of 15-20 minutes each, combined with light acceleration. Using the RS gel alone enables patients to overcome any hypersensitivity in a very short time after treatment and does not require the use of the accelerator. The result is excellent whitening results coupled with little to no sensitivity for those undergoing the treatment. **FOR BEST RESULTS, USE WITH BEYOND® ACCELERATOR.**

14



## BEYOND® MAX5® Solution Treatment Kit

- Features our popular 35% H<sub>2</sub>O<sub>2</sub> W Gel 1 formula whitening gel
- Our specially formulated RS Gel contains special ingredients to reduce the risk of sensitivity
- The combination of the two gels rebuilds, strengthens, and whitens tooth enamel
- Gel's proprietary formula improves oral health while whitening
- Produces dramatic glossy, white results
- Clinically proven results

### Package Contents

- 5 Dual-barrel syringes of Max5 Solution W gel (2.6g/2.2 mL)
- 5 Mixing tips of Max5 Solution W gel syringe
- 5 Pouches of Max5 Solution RS gel (11g/9.5mL)
- 5 Spatulas
- 5 Brushes
- 5 Mixing jars
- 5 Face protection cloths
- 5 Vita shade guides
- 5 After-treatment guides
- 1 Instruction manual



## BEYOND® SOLO® Solution Treatment Kit

- Features our popular 35% H<sub>2</sub>O<sub>2</sub> W Gel 1 formula whitening gel
- Our specially formulated RS Gel contains special ingredients to reduce the risk of sensitivity
- Gel's proprietary formula improves oral health while whitening
- Produces dramatic glossy, white results
- Compact packaging, economical and practical
- Individualized packaging for single patients



BEYOND® Solo® Solution was created with patients with known cases of sensitive teeth in mind. Solo® Solution consists of a W Gel of 35% H<sub>2</sub>O<sub>2</sub> dual barrel syringe, which is then mixed with a RS Gel that contains known desensitizing ingredients such as potassium nitrate and restoring ingredients such as hydroxyapatite. This reduces the concentration of hydrogen peroxide to below 6% and helps aid in reduction of sensitivity as well as remineralization of the enamel. The result is excellent whitening results coupled with little to no sensitivity for those undergoing the treatment. **FOR BEST RESULTS, USE WITH BEYOND® ACCELERATOR.**

## Our Most Economical, Single Patient Whitening Kit

### Package Content

- 1 Dual-barrel syringe of Solo Solution W gel (2.6g/2.2 mL)
- 1 Mixing tip for Solo Solution W gel syringe
- 1 Pouch of Solo Solution RS gel (11g/9.5mL)
- 1 Face protection cloth
- 1 Spatula
- 1 Brush
- 1 Vita shade guide
- 1 Instruction manual
- 1 After-treatment guide





## BEYOND® MAGNA™ Solution Treatment Kit

- Features our popular 35% H<sub>2</sub>O<sub>2</sub> W Gel 1 formula whitening gel
- Our specially formulated RS Gel contains special ingredients to reduce the risk of sensitivity
- Gel's proprietary formula improves oral health while whitening
- Produces dramatic glossy, white results
- Contains special ingredients to reduce the risk of sensitivity
- Clinically proven results
- Compact packaging, a full size two patient treatment kit

BEYOND® Magna™ Solution was created with patients with known cases of sensitive teeth in mind. Magna™ Solution consists of a W Gel of 35% H<sub>2</sub>O<sub>2</sub> dual barrel syringe, which is then mixed with a RS Gel that contains known desensitizing ingredients such as potassium nitrate and restoring ingredients such as hydroxyapatite. This reduces the concentration of hydrogen peroxide to below 6% and helps aid in reduction of sensitivity as well as remineralization of the enamel. The result is excellent whitening results coupled with little to no sensitivity for those undergoing the treatment. **FOR BEST RESULTS, USE WITH BEYOND® ACCELERATOR.**

### Package Contents

- 2 Magna Solution W gel 1 dual-barrel syringes (4.0g/3.5 mL)
- 2 Mixing tips of Magna Solution W gel dual-barrel syringe
- 2 Pouches of Magna Solution RS gel (16g/13.8mL)
- 2 Brushes
- 2 Spatulas
- 2 Face protection cloths
- 2 Vita shade guides
- 1 Instruction manual

## Economical Full-Size, Twin-Pack Whitening Kit



## Package Content

- 1 Jar of Blast® Solution silicon dioxide catalyst powder (0.7g)
- 1 Bottle of Blast® Solution 35% hydrogen peroxide aqueous solution (2mL)
- 1 Pouch of Blast® Solution RS gel (16g/13.8mL)
- 1 Spatula
- 1 Brush
- 1 Face protection cloth
- 1 Vita shade guide
- 1 Instruction manual

# Special Designed Formula For Tetracycline and Severe Stains

BEYOND® Blast® Solution was created with patients with known cases of severe or tetracycline stains in mind. Blast® Solution consists of a W Gel of 35% H<sub>2</sub>O<sub>2</sub> aqueous solution, which is mixed with silicon dioxide catalyst powder into a wet paste. The paste is then mixed with a RS Gel that contains known desensitizing ingredients such as potassium nitrate and restoring ingredient such as hydroxyapatite. This reduces the concentration of hydrogen peroxide to below 6% and helps aid in reduction of sensitivity as well as remineralization of the enamel. The result is excellent whitening results coupled with little to no sensitivity for those undergoing the treatment. **FOR BEST RESULTS, USE WITH BEYOND® ACCELERATOR.**

## BEYOND® Blast® Solution Treatment Kit

- 35% hydrogen peroxide manual-mix whitening paste
- Our specially formulated RS Gel contains special ingredients to reduce the risk of sensitivity
- Formulated to reduce the appearance of difficult and severe stains
- Effective pre-treatment whitening boost for tetracycline and fluoride stains
- For use with any BEYOND® Whitening Kit or light system
- Individualized packaging for single patients
- Clinically proven results







## 6% H<sub>2</sub>O<sub>2</sub> Premium Whitening Formula

This 6% Hydrogen Peroxide (H<sub>2</sub>O<sub>2</sub>) Teeth Whitening Kit is designed to deliver professional results in just 36 minutes. The OSMO Five -Patient Whitening Kit features a safe, advanced formula gel that can be used in-office without gingival protection. This kit includes material for five patients. A typical treatment is three consecutive cycles of 12 minutes each, designed to work in conjunction with light acceleration. The dual barrel 6% syringe in this kit includes enough teeth bleaching gel for three chairside cycles.

### Package Contents

- 5 Dual-barrel syringes of Osmo whitening gel (2.6g/2.2 mL)
- 5 Mixing tips of Osmo gel dual-barrel syringe
- 5 Face protection cloths
- 5 Vita shade guides
- 5 After-treatment guides
- 1 Instruction manual

### BEYOND<sup>®</sup> OSMO<sup>®</sup> Treatment Kit

- Specially formulated 6% H<sub>2</sub>O<sub>2</sub> gel designed for excellent accelerated results
- Gel's gentle formula does not require gingival protection
- Static mixing tip provides a one-step gel application
- Brilliant results in only 36 minutes
- Guaranteed to deliver shiny, white results with NO sensitivity



## BEYOND® OSMO® Teeth Whitening Pen Twin Pack

- Provides a renewed white smile every day
- Safe for daily use and whitening maintenance
- Gentle brush-on pen for easy application
- Can be used alone or in conjunction with all BEYOND® whitening procedures
- Enough material for 30 days
- Mint flavored 6% H<sub>2</sub>O<sub>2</sub> gel for a white and fresh smile



This pack features two teeth whitening pens, each with enough hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) gel to last for 30 days. Use the pen for just minutes daily as maintenance following a professional teeth whitening service or as a stand alone at-home teeth whitening product.

### Package Content

- 2 Osmo® whitening pens (3.2g/2.8 mL)

# Premium Formula for Home Whitening Maintenance



# Revolutionize your teeth in just 1 Hour

## BEYOND® Gemini® Teeth Whitening Accelerator Home Edition

- Next generation LED for maximum acceleration
- No sensitivity
- Removes years of stains caused by coffee, red wine and smoking. 6 to 8 shade improvements
- Embedded magnets line light to tray for a perfect fit every time
- Simply rechargeable via USB or outlet
- Clinically proven to kill germs
- Refill kits available for continued maintenance



The Gemini® home whitening system guarantees fresh, safe and brilliant results in just one hour. The gentle 6% hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) formula is a perfect at home teeth bleaching solution. With over 20 years of experience, BEYOND® has generated an advanced system that is designed for you and defined for your smile in the comfort of your home. Gemini® Teeth Whitening Accelerator includes whitening gel for five – twenty minute sessions, Gemini® LED light unit, Gemini tray, rechargeable lithium battery, charge cord and USB adaptor and bleaching shade guide.

### Package Contents

- 1 Gemini® LED light unit
- 5 Syringes of Gemini® whitening gel (1.9g/1.7mL)
- 5 Tips of Gemini® whitening gel syringe
- 1 Gemini® tray
- 1 Gemini® power unit with rechargeable lithium battery
- 1 Gemini® charge cord
- 1 Gemini® USB/electrical adaptor
- 1 Bleaching shade guide
- 1 Instruction manual



## Package Content

- 6 syringes CoreWhite gel (2.36g/2.2 mL)
- 2 Sheets of EVA impression material
- 1 CoreWhite tray box
- 1 Manual booklet
- 1 Carrying case

# Take-home tray whitening available through the dental office

## BEYOND® CoreWhite® Take-Home Whitening Kit

- 30 minutes a day for 14 days – no messy overnight whitening
- Visibly whiter teeth in five days
- Gentle hydrogen peroxide formula minimizes risk of sensitivity (6% H<sub>2</sub>O<sub>2</sub> gel)
- EVA material included for a comfortable customized tray
- Mint and Orange flavors
- CoreWhite® Refill Kits also available



New and gentle 6% hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) formula guarantees a beautiful white smile. Offered through dental offices, the CoreWhite take-home tray teeth whitening system provides professional results from the comfort of your home for only 30 minutes a day. This kit includes material to make custom trays, and a handy carrying case with bag. A great option to have your teeth whitened in a simple, quick, and inexpensive manner.



# Advanced, Daily Whitening Maintenance



22

## Pearl White® Whitening Toothpaste

- Unique formula includes calcium peroxide that acts as an additional whitening agent
- Provides long lasting fresh breath protection
- Advanced tartar control
- No animal testing, no animal ingredients
- Made in the USA

Developed by a cosmetic dentist, Pearl White® toothpaste is designed to both gently brush away surface stains as well as provide tartar control. Our formula contains micro-crystals that burst and whiten teeth while you brush. Our teeth whitening toothpaste contains quality ingredients freshens breath while supporting enamel mineralization that enables superior whitening results.

Available in 3 unique formulas: **Advanced, Sensitivity, Green Tea** and two sizes: **Standard 4.2oz / 120g, Travel 1.4oz / 40g**





NET WT 4.2 OZ (120ml)

PEARL WHITE

Whitening Toothpaste  
Sensitivity Formula