

Bombillas para Equipos Dentales.

Fabricant	Modèle	Réf. °K	Tensión(V)	potencia(W)	Composition	DegreK	Tecnomed
Adec	6300 (300 h)	°K203	24	150		11407	88820
	6300 (40 h)	°K209	24	150		11410	88860
	Performer	°K208	12	55		11409	88850
Anthos	Venus	°K233	17	95		11415	97812
Belmont	–	°K208	12	55		11409	88850
Castellini	Luna	°K217	12	75		11411	88910
	Monta	°K212	12	50		88890	
	Solar	°K212	12	50		88890	
	Vega	°K212	12	50		88890	
Eurodent	Isolite (300 h)	°K203	24	150		11407	88820
	Isolite (40 h)	°K209	24	150		11410	88860
	Ruinito (300 h)	°K203	24	150		11407	88820
	Ruinito (40 h)	°K209	24	150		11410	88860
Faro	70 (300 h)	°K203	24	150		11407	88820
	70 (40 h)	°K209	24	150		11410	88860
	B4843 (300 h)	°K203	24	150		11407	88820
	B4843 (40 h)	°K209	24	150		11410	88860
	EDI	°K233	17	95		11415	97812
	Forest	°K234	24	150	avec culot ceramique	88950	
	Koolbeam	°K234	24	150	avec culot ceramique	88950	
	S2000 (tranfo. FARO 17V)	°K204	17	95	avec culot ceramique	11408	88830
	S70 (tranfo. FARO de 17V)	°K204	17	95	avec culot ceramique	11408	88830
	S90 (tranfo. FARO de 17V)	°K204	17	95	avec culot ceramique	11408	88830
	Sunlight (300 h)	°K203	24	150		11407	88820
	Sunlight (40 h)	°K209	24	150		11410	88860
Fedesa	Midway	°K204	17	95	avec culot ceramique	11408	88830
	Scorpio	°K204	17	95	avec culot ceramique	11408	88830
	lineo	°K204	17	95	avec culot ceramique	11408	88830
Gallus	Alpha 1975	°K201	0	200	6,6A		
	Delta	°K202	0	200	6,6A	11406	
	G308	°K202	0	200	6,6A	11406	
Healthco	Aristocrat	°K208	12	55		11409	88850
	Celebrity	°K208	12	55		11409	88850
Kavo	1400 (culot 2 broches, 300h)	°K203	24	150		11407	88820
	1400 (culot 2 broches, 40h)	°K209	24	150		11410	88860
Philips	Oralux	°K212	12	50		88890	
Planmeca	Delight (300 h)	°K203	24	150		11407	88820
	Delight (40 h)	°K209	24	150		11410	88860
Quetin	–	°K208	12	55		11409	88850
	White light	°K201	0	200	6,6A		

Bombillas para Equipos Dentales.

Fabricant	Modèle	Réf. °K	Tensión(V)	potencia(W)	Composition	DegreK	Tecnomed
Ritter	D428 (300 h)	°K203	24	150		11407	88820
	D428 (40 h)	°K209	24	150		11410	88860
	D44 (culot 2 broches, 300 h)	°K203	24	150		11407	88820
	D44 (culot 2 broches, 40h)	°K209	24	150		11410	88860
	Lumostar new (300 h)	°K203	24	150		11407	88820
	Lumostar new (40 h)	°K209	24	150		11410	88860
	Ritter 130 (300 h)	°K203	24	150		11407	88820
	Ritter 130 (40 h)	°K209	24	150		11410	88860
Sirona	Fantastic (n° série de 0 à 29999)	°K220	24	150		88920	
	Sirolux	°K208	12	55		11409	88850
Villa	Moonlight	°K208	12	55		11409	88850
Vitali	Vitali Cammbologna	°K204	17	95	avec culot ceramique	11408	88830
autres	Eurocad	°K224	13.8	50		11413	
	Eurostar	°K224	13.8	50		11413	
	Eurostar	°K230	12	50		11414	88940
	Magnat X,	°K212	12	50		88890	
	Maruchi Eight	°K212	12	50		88890	
	Virginia (Chayes)	°K208	12	55		11409	88850