

# CEFLA DENTAL HOSES

The dental hoses for Tecnomed Italia dental units are made of silicone; they are Latex Free and comply with the UNI CEI EN ISO 13485 quality system.

The chemically stable gold-plated contacts ensure high electrical conductivity and maximum handpiece performance.

The supply tubes comply with REGULATION (EU) No. 10/2011 concerning plastic materials objects intended to meet food products and the exhaust tubes are resistant to the lubrication oils.

Tecnomed Italia dental hoses meet the strictest medical standards, complying with the essential requirements of the European Directives (93/42 / EEC) 2007/47 / EC on medical devices and 2000/95 / EC RoHs on the restriction of the use of certain specific substances in electrical and electronic equipment.

A complete range of dental hoses:

- Turbines
- Micromotors
- Scalers
- Syringes
- Freely Adaptable Dental Hoses










## Features

- ⊕ High quality silicone sheath.
- ⊕ Latex free.
- ⊕ 12 months warranty










## Technical data

- ⊕ Maximum pressure supply/return air 4,5 bar / 65 psi
- ⊕ Maximum pressure air/water spray 4 bar / 58 psi
- ⊕ Tensile strenght 70N










## MICROMOTORS TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MICROMOTOR MC2	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9  STERN S200, 220, 250, 280, 300, 320, 380  CASTELLINI SKEMA 5, 6, 8 PUMA ELI				PC707CF
MICROMOTOR MC2	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9  STERN S200, 220, 250, 280, 300, 320, 380  CASTELLINI SKEMA 5, 6, 8 PUMA ELI				PC709CF
MICROMOTOR MC3 M40	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9  STERN S200, 220, 250, 280, 300, 320, 380  CASTELLINI SKEMA 5, 6, 8 PUMA ELI				PC710CF




# MICROMOTORS TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MICROMOTOR  MC3 M40	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9  STERN S200, 220, 250, 280, 300, 320, 380  CASTELLINI SKEMA 5, 6, 8 PUMA ELI		 <p>Spray regulation Swivel Safety Handpiece</p> <p>ITAL GRAY RAL7038</p> <p>1160 mm - Ø 10,5 mm</p> <p>1250 mm</p>		PC710CFSH
MICROMOTOR  MC3 M40	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9  STERN S200, 220, 250, 280, 300, 320, 380  CASTELLINI SKEMA 5, 6, 8 PUMA ELI		 <p>Spray regulation Swivel</p> <p>ITAL GRAY RAL7038</p> <p>1350 mm - Ø 10,5 mm</p> <p>1440 mm</p>		PC712CF
MICROMOTOR  MC3 M40	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9  STERN S200, 220, 250, 280, 300, 320, 380  CASTELLINI SKEMA 5, 6, 8 PUMA ELI		 <p>Spray regulation Swivel Safety Handpiece</p> <p>ITAL GRAY RAL7038</p> <p>1350 mm - Ø 10,5 mm</p> <p>1440 mm</p>		PC712CFSH

# MICROMOTORS TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MICROMOTOR  iXR3L iMMrL	ANTHOS A3, 5/+, 6/+, 7+, 9, L6, 9, R7  STERN S200, 250, 280TRc, 380TRc, 300, 220TR, 230TR  CASTELLINI SKEMA 5, 6, 8				PC733CX3
MICROMOTOR  iXR3L iMMrL	ANTHOS A3, 5/+, 6/+, 7+, 9, L6, 9, R7  STERN S200, 250, 280TRc, 380TRc, 300, 220TR, 230TR  CASTELLINI SKEMA 5, 6, 8				PC733CX3L
MICROMOTOR  iXS4 iMMs	ANTHOS A7+ L9, R7  STERN S320TR, 280TRC, 300, 380TRC  CASTELLINI SKEMA 5, 6, 8				PC744CX4 BEFORE 2016  PC744CX4N SINCE 2016









# MICROMOTORS TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGTH / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MICROMOTOR  iXS4 iMMs	ANTHOS A7+ L9, R7  STERN S320TR, 280TRC, 300, 380TRC  CASTELLINI SKEMA 5, 6, 8		 <p>Spray regulation</p> <p>Swivel</p> <p>ITAL GRAY RAL7038</p> <p>1350 mm - Ø 10 mm</p> <p>1440 mm</p>		PC744CX4L BEFORE 2016  PC744CX4LN SINCE 2016










# HIGHSPEED TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGTH / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
HIGHSPEED HSH	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9  STERN S200, 220, 250, 280, 300, 320, 380  CASTELLINI SKEMA 5, 6, 8 PUMA ELI		<p>1160 mm - Ø 10,5 mm 1250 mm</p>		PC704CF
HIGHSPEED HSH	SKEMA 4 DUO		<p>1160 mm - Ø 10,5 mm 1250 mm</p>		PC704CFW
HIGHSPEED HSH	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9  STERN S200, 220, 250, 280, 300, 320, 380  CASTELLINI SKEMA 5, 6, 8 PUMA ELI		<p>1350 mm - Ø 10,5 mm 1440 mm</p>		PC706CF










# HIGHSPEED TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
HIGHSPEED HSH	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9				PC704CFSH
	STERN S200, 220, 250, 280, 300, 320, 380  CASTELLINI SKEMA 5, 6, 8 PUMA ELI				
HIGHSPEED HSH	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9				PC706CFSH
	STERN S200, 220, 250, 280, 300, 320, 380  CASTELLINI SKEMA 5, 6, 8 PUMA ELI				










# SYRINGE TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
<p>SYRINGE</p> <p>6F LUZZANI O.F. MINIBRIGHT INOX</p>	<p>ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9</p> <p>STERN S200, 220, 250, 280, 300, 320, 380</p> <p>CASTELLINI SKEMA 5, 6, 8 PUMA ELI</p>		<p>ITAL GRAY RAL7039</p> <p>1160 mm - Ø 10,5 mm</p> <p>1250 mm</p> 		<p>PC714EL</p>
<p>SYRINGE</p> <p>3F LUZZANI MINILIGHT INOX</p>	<p>ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9</p> <p>STERN S200, 220, 250, 280, 300, 320, 380</p> <p>CASTELLINI SKEMA 5, 6, 8 PUMA ELI</p>		<p>ITAL GRAY RAL7038</p> <p>1160 mm - Ø 10,5 mm</p> <p>1250 mm</p> 		<p>PC714C</p>
<p>SYRINGE</p> <p>3F LUZZANI MINILIGHT INOX</p>	<p>ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9</p> <p>STERN S200, 220, 250, 280, 300, 320, 380</p> <p>CASTELLINI SKEMA 5, 6, 8 PUMA ELI</p>		<p>ITAL GRAY RAL7038</p> <p>1350 mm - Ø 10,5 mm</p> <p>1440 mm</p> 		<p>PC714L</p>










# SCALER TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
SCALER EMS	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9  STERN S200, 220, 250, 280, 300, 320, 380  CASTELLINI SKEMA 5, 6, 8 PUMA ELI				PC701CF
SCALER EMS	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9  STERN S200, 220, 250, 280, 300, 320, 380  CASTELLINI SKEMA 5, 6, 8 PUMA ELI				PC703CF
SCALER PZ5	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9  STERN S200, 220, 250, 280, 300, 320, 380  CASTELLINI SKEMA 5, 6, 8 PUMA ELI				PC701CP5

# SCALER TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
SCALER PZ5	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9  STERN S200, 220, 250, 280, 300, 320, 380  CASTELLINI SKEMA 5, 6, 8 PUMA ELI				PC701CP5F
SCALER EMS O. F.	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9  STERN S200, 220, 250, 280, 300, 320, 380  CASTELLINI SKEMA 5, 6, 8 PUMA ELI				PC701CFE
SCALER EMS O. F.	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9  STERN S200, 220, 250, 280, 300, 320, 380  CASTELLINI SKEMA 5, 6, 8 PUMA ELI				PC703CFEL










# SCALER TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
SCALER PIEZOLIGHT 6 OLD PZ7-SCA3	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI		 ITAL GRAY RAL7038 1160 mm - Ø 8 mm 1250 mm		PC701CF4
SCALER PIEZOLIGHT 6 OLD PZ7-SCA3	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI		 ITAL GRAY RAL7038 1350 mm - Ø 8 mm 1440 mm		PC703CF4
SCALER SATELEC	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI		 ITAL GRAY RAL7038 1160 mm - Ø 7 mm 1250 mm		PC701CFS







# SCALER TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
SCALER SATELEC	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9  STERN S200, 220, 250, 280, 300, 320, 380  CASTELLINI SKEMA 5, 6, 8 PUMA ELI		<p>Spray regulation Swivel</p> <p>ITAL GRAY RAL7038</p> <p>1350 mm - Ø 7 mm</p> <p>1440 mm</p>		PC703CFS
SCALER SATELEC NEWTRON LED	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9  STERN S200, 220, 250, 280, 300, 320, 380  CASTELLINI SKEMA 5, 6, 8 PUMA ELI		<p>Spray regulation Swivel</p> <p>ITAL GRAY RAL7038</p> <p>1160 mm - Ø 8 mm</p> <p>1250 mm</p>		PC701CFSL
SCALER SATELEC NEWTRON LED	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9  STERN S200, 220, 250, 280, 300, 320, 380  CASTELLINI SKEMA 5, 6, 8 PUMA ELI		<p>Spray regulation Swivel</p> <p>ITAL GRAY RAL7038</p> <p>1350 mm - Ø 8 mm</p> <p>1440 mm</p>		PC703CFSL

# SCALER TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
SCALER PIEZOSTERIL 5 SURGISON	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9  STERN S200, 220, 250, 280, 300, 320, 380  CASTELLINI SKEMA 5, 6, 8 PUMA ELI				PC705CPS
SCALER PIEZOLIGHT 6 PZ7-SCA3 NEW	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9  STERN S200, 220, 250, 280, 300, 320, 380  CASTELLINI SKEMA 5, 6, 8 PUMA ELI				PC741C
SCALER PIEZOLIGHT 6 PZ7-SCA3 NEW	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9  STERN S200, 220, 250, 280, 300, 320, 380  CASTELLINI SKEMA 5, 6, 8 PUMA ELI				PC741CL

# SCALER TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
SCALER  PIEZOLIGHT 6 PZ7-SCA3 NEW	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9  STERN S200, 220, 250, 280, 300, 320, 380  CASTELLINI SKEMA 5, 6, 8 PUMA ELI				PC742C
SCALER  PIEZOLIGHT 6 PZ7-SCA3 NEW	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9  STERN S200, 220, 250, 280, 300, 320, 380  CASTELLINI SKEMA 5, 6, 8 PUMA ELI				PC742CL