

FICHA TÉCNICA

GUANTES DE LÁTEX CON POLVO BESTDENT		
REFERENCIAS	24684: Tamaño XL 24685: Tamaño L 24686: Tamaño M 24687: Tamaño S 24688: Tamaño XS	
FABRICANTE / RESPONSABLE	TOP GLOVE SDN. BHD. Lot 4969, Jalan Teratai, Batu 6, Off Jalan Meru 41050 Klang, Selangor D.E., Malaysia	
CLASIFICACIÓN	Producto Sanitario Clase I, Reglamento 2017/745/UE. EPI Categoría III, según el Reglamento 2016/425/UE.	
DESCRIPCIÓN	Guantes de látex con polvo, no estériles, de un solo uso. Ambidiestros.	
COMPOSICIÓN	100% goma natural.	
ESPECIFICACIONES	Tamaños: L, M, S y XS	
EMBALAJE	Caja con 100 unidades	
CADUCIDAD	5 años.	
FINALIDAD PREVISTA	Protección de barrera de la mano del usuario para la realización de exámenes médicos o clínicos, procedimientos diagnósticos y terapéuticos, en odontología y también para fines de laboratorio.	
INFORMACIÓN SOBRE LA SEGURIDAD	Los guantes están fabricados con látex, que puede causar reacciones alérgicas en algunas personas. Si esto ocurre, discontinuar su uso inmediatamente y consultar con el médico.	
ALMACENAMIENTO Y TRANSPORTE	Almacenar en lugar fresco y seco. Mantener alejado de la luz solar. Proteger de la humedad y evitar el calor excesivo.	

TECHNICAL DATA SHEET

POWDERED LATEX GLOVES BESTDENT		
REFERENCES	24684: Size XL 24685: Size L 24686: Size M 24687: Size S 24688: Size XS	
MANUFACTURER	TOP GLOVE SDN. BHD. Lot 4969, Jalan Teratai, Batu 6, Off Jalan Meru 41050 Klang, Selangor D.E., Malaysia	
CLASSIFICATION	Class I Medical Device, Regulation 2017/745/EU. PPE Category III, according to Regulation 2016/425/EU.	
DESCRIPTION	Powdered latex gloves, non-sterile, single use. Ambidextrous.	
COMPOSITION	100% natural rubber.	
SPECIFICATION	Sizes: L, M, S and XS	
PACKAGE	Box with 100 units.	
SHELF LIFE	5 years.	
INTENDED USE	Barrier protection of the user's hand for medical or clinical examinations, diagnostic and therapeutic procedures, in dentistry and also for laboratory purposes.	
SAFETY INFORMATION	The gloves are made of latex, which may cause allergic reactions in some people. If this occurs, discontinue use immediately and consult a physician.	
STORAGE AND TRANSPORT	Store in a cool and dry place. Keep away from sunlight. Protect from moisture and avoid excessive heat.	